



Zumax Medical
Beyond Vision

Beyond
your visual limits!



CE FDA



THROUGHT-THE-LENS (TTL) Loupes



- 2.5 and 3.0 magnifications available
Extremely high resolution
- Precision Optics
- Custom Fit and Custom Working Distance
- Greatest Brightness and Depth of Field
- Ultra-lightweight, as little as 2 oz.(57g) (including the frame).



All lenses are multi-coated and anti-reflective, combining with infinity corrected optical system.



Soft silicone nose pad.



Two kinds of frames for your option which provides the best comfort at work.
1. Sport Frame, three colors for your optional.
2. Aluminium magnesium alloy Frame, black



All loupes include a protective hard case, microfiber wipe cloth, neck strap.

Specification :

Magnification	Working distance	Field of view
2.5x	340mm/13"	60mm/2.4"
2.5x	420mm/17"	70mm/2.8"
2.5x	460mm/18"	80mm/3.1"
2.5x	500mm/20"	95mm/3.8"

Magnification	Working distance	Field of view
3.0x	340mm/13"	40mm/1.6"
3.0x	420mm/17"	46mm/1.8"
3.0x	460mm/18"	50mm/2.0"
3.0x	500mm/20"	54mm/2.1"

ZUMAX MEDICAL CO., LTD.

36 Zhiying Street, Suzhou New District
Jiangsu, China 215129
Tel: +86 512 66650100 ext. 502/504/510
Fax: +86 512 66909655
Email: sales@zumaxmedical.cn
www.zumaxmedical.com