

ROTARY
GRINDING &
POLISHING
INSTRUMENTS

VRTEČI BRUSILNI & POLIRNI INSTRUMENTI



MADE IN GERMANY

DENTISTRY

PAGE STRAN		ABRASIVES ABRAZIVI	3D Printing Hybrid Material Hibridni Materiali za 3D Tiskanje						Gold / Amalgam Zlato / Amalgam		NEM Alloys NEM Kovine		Titanium Titan		Acrylics Akriat		Prophylaxis Profilaksa		PROCESSING STEPS STOPENJSKI PROCES	GRIT GROBOST					
CERASTONE green	8	SIC																						F	
CERASTONE white	8	ALOX																						EF	
DIACOMP® plus	9																			2			M	F	
EASYCOMP	10	SIC																		2			M	F	
ECOCOMP	11	SIC																		2			M	F	
COMPOSOF	11	ALOX																		3		C	M	F	
DIAPOL®	12																			3		C	M	F	
DIACERA	13																			2			M	F	
DIAPRO	14																			2			M	F	
ECOCERAM	15	SIC																		3		C	M	F	
UNIVERSAL	15	SIC																		1			M		
EVEFLEX	16	SIC																		3		C	M	F	
FLEXI-DIA	17																			4		C	M	F	EF
FLEXI-D	17	ALOX																		4		C	M	F	EF
EVEBRUSH	18	SIC																		1				F	
BRACKET plus	18	SIO																							
PROPHYCLEAN	19	SIO																		2		C	M		
DIASTRIPS	19																			3		C	M	F	

particularly suitable
zelo primerni

suitable
primerni

VC very coarse
zelo grobi

C coarse
grobi

M medium
srednji

F fine
fini

EF extra fine
ekstra fini

Diamond impregnated
Impregnirani z diamantom

SIC Silicon carbide
Silicijev karbid

ALOX Aluminium oxide
Aluminijev oksid

SIO Silicon oxide
Silicijev oksid

Disposable
Za enkratno uporabo

Recommended RPM
Priporočljivi vrtljaji

Packaging unit
Enota pakiranja

TWIST

ONE SHAPE

FOR EVERY SURFACE

ENA OBLIKA ZA VSE POVRŠINE



**One instrument
for every dental surface**
En polirnik za vse površine zoba



**Preservation of the
surface structure**
Ohranjanje površinske strukture



Extraordinary durability
Izjemna vzdržljivost



Reduced heat generation
Manj segrevanja




Intuitive application
Intuitivna uporaba



Brilliant polishing results
Brilijantni rezultati poliranja

THE OCCLUSION SOLUTION

OKLUZIJSKA REŠITEV



Designed specifically for occlusal polishing
Dizajniran posebej za okluzijska poliranja



Applicable without polishing paste
Uporaba brez polirne paste



Intuitive application
Intuitivna uporaba




Reaches deep fissures due to twisted polishing bristles
Zvite polirne konice dosežejo globoke fisure



Easily reprocessable brush shape
Enostavno dodatno procesiranje



Brilliant polishing results
Brilijantni rezultati poliranja



Why I am working with TWIST & OCCLUFLEX

The surface quality of my direct composite restorations, posterior and anterior, changed radically since I started using EVE DIACOMP *plus*. In only two polishing steps, it is possible to get a remarkably smooth and shiny surface result that not only enhances the aesthetics, but also significantly improves the longevity of my direct composite restorations.

The application of the TWIST and OCCLUFLEX shapes just differs in the area of application in the oral cavity. TWIST polishers are more adequate for anterior restorations with a flatter surface. On the other hand, OCCLUFLEX is more suitable for posterior restorations due to its design. The special geometry of the latter ones make them more efficient when polishing in deeper regions of my restorations.

Dr. Vasco Nunes da Silva
Dentist and lecturer

Zakaj uporabljam TWIST & OCCLUFLEX

Odkar sem začel uporabljati EVE DIACOMP *plus*, se je kakovost površin pri mojih direktnih anteriornih in posteriornih kompozitnih restavracijah korenito spremenila. S samo dvema korakoma poliranja je sedaj mogoče doseči izredno gladko in sijočo površino, ki ne le izboljša estetiko, ampak tudi bistveno podaljša življensko dobo restavracij.

Oblika TWIST in OCCLUFLEX polirnikov se razlikuje glede na področje obdelave pri uporabi v ustni votlini. Polirniki TWIST so primernejši za anteriorne restavracije, z bolj ravno površino. Zaradi svoje oblike, pa so OCCLUFLEX polirniki primernejši za posteriorne restavracije. Med poliranjem globokih predelov restavracij, naredi posebna geometrija polirnih konic polirnike OCCLUFLEX še učinkovitejše.

Dr. Vasco Nunes da Silva
Zobozdravnik in predavatelj



Dr. Vasco Nunes da Silva
Dentist and lecturer
Zobozdravnik in predavatelj



DIACOMP® plus | TWIST

Page 9 / Stran 9



DIACOMP® plus | OCCLUFLEX

Page 9 / Stran 9




CERASTONE

Composite, Feldspar Ceramic, Amalgam, Gold
Kompozit, Leucitna Keramika, Amalgam, Zlato

GRINDING SYSTEM | BRUSILNI SISTEM

- Ceramic bonded abrasives
- Silicon carbide and corundum grit
- For detailed shaping of dental materials
- Abrazivi z vezano keramiko
- Silicijev karbid in korundna zrnca
- Za detajlno oblikovanje zobozdravstvenih materialov



Application | Uporaba:

Grit | grobost:

Bonding | vezava:

Silicon carbide abrasives for Feldspar Ceramic, Amalgam and Gold Silicijevi karbidini polirniki za Leucitno Keramiko, Amalgam in Zlato	FINE fini	MEDIUM srednji	RA
Arkansas stones for Composite finishing Arkansas kamni za končno obdelavo kompozita	EXTRA FINE ekstra fini	HARD trdi	RA
Arkansas stones for Composite finishing Arkansas kamni za končno obdelavo kompozita	EXTRA FINE ekstra fini	HARD trdi	FG

FINE FINI



15.000 - 25.000 min⁻¹
max. 30.000 min⁻¹



10, 100 pieces / kosov



Item / Artikel:	CSG-RA-601	CSG-RA-645	CSG-RA-649	CSG-RA-661	RA
REF:	8601	8603	8604	8605	
Dimensions / Dimenzija (mm):	3 x 3	3 x 7	2,5 x 6	2,5 x 7	

EXTRA FINE EKSTRA FINI



RA: 15.000 - 25.000 min⁻¹
max. 30.000 min⁻¹



FG: 80.000 - 100.000 min⁻¹
10, 100 pieces / kosov



Item / Artikel:	CSW-RA-601	CSW-RA-645	CSW-RA-649	CSW-RA-661	RA
REF:	8631	8633	8634	8636	
Dimensions / Dimenzija (mm):	CSW-FG-601	CSW-FG-645	CSW-FG-649	CSW-FG-661	FG
REF:	8694	8696	8697	8698	
Dimensions / Dimenzija (mm):	3 x 3	3 x 7	2,5 x 6	2,5 x 7	

DIACOMP® plus

Composite, Hybrid Ceramic, 3D Printing Hybrid Material
Kompozit, Hibridna Keramika, Hibridni Materiali za 3D Tiskanje



DIAMOND IMPREGNATED POLISHING SYSTEM
Z DIAMANTI IMPREGNIRAN POLIRNI SISTEM

- 2-step polishing system for high-shine on all composite materials
- Polishers specifically adapted to each processing step
- Best results even on very hard composites
- 2-stopenjski polirni sistem za visok sijaj na vseh kompozitih
- Polirniki prilagojeni posamezni stopnji procesa
- Najboljši rezultati tudi na zelo trdih kompozitih

Set RA 345 REF 9095

DCP-Ofm | DCP-Off



Set RA 342 REF 9086

DT-DCP14m | DT-DCP14f



Set RA 344 REF 9094

DCP-W11m | DCP-Ofm | DCP3m | DCP2m
DCP-W11f | DCP-Off | DCP3f | DCP2f



MEDIUM Pre-polishing
SREDNJI Predpoliranje



FINE High-shine polishing
FINI Poliranje visoki sijaj



2 PROCESSING STEPS
STOPENJSKI PROCES



3.000 - 8.000 min⁻¹
max. 20.000 min⁻¹



1, 10 piece / kosov

RA	DT-DCP14m	DCP-Ofm	DCP-W11m	DCP1m	DCP3m	Item / Artikel:
	6184	6175	6181	6171	6173	REF:
	14 x 1,6	5,5 x 8,8	3 x 7,5	4 x 10	6 x 7,5	Dimensions / Dimenzija (mm):



3.000 - 8.000 min⁻¹
max. 20.000 min⁻¹



1, 10 piece / kosov

RA	DT-DCP14f	DCP-Off	DCP-W11f	DCP1f	DCP3f	Item / Artikel:
	6284	6275	6281	6271	6273	REF:
	14 x 1,6	5,5 x 8,8	3 x 7,5	4 x 10	6 x 7,5	Dimensions / Dimenzija (mm):

EASYCOMP

Composite · Kompozit

POLISHING SYSTEM | POLIRNI SISTEM



- Economic single-use polishing system for composite
- Elimination of reprocessing
- Small geometries simplify the accessibility of difficult to reach areas
- Optimised range of available shapes for all surface structures
- Ekonomičen polirni sistem za enkratno uporabo na kompozitih
- Ni potrebno dodatno procesiranje
- Manjše geometrije omogočajo enostavnejši dostop do težje dostopnih predelov
- Optimiziran nabor razpoložljivih oblik za vse površinske strukture

2 PROCESSING STEPS STOPENJSKI PROCES



MEDIUM Pre-polishing
SREDNJI Predpoliranje



FINE Polishing
FINI Poliranje



3.000 - 8.000 min⁻¹
* 3.000 - 5.000 min⁻¹
max. 20.000 min⁻¹



40 pieces / kosov

Set EY-10

20x EY-10m | 20x EY-10f



Set EY-OF *

20x EY-OFm | 20x EY-OFf



Item / Artikel:	Set EY-10	EY-10m	EY-10f	Set EY-OF	EY-OFm	EY-OFf	RA
REF:	9033	3280	3380	9032	3275	3375	
Dimensions / Dimenzija (mm):	—	11 x 1,6	11 x 1,6	—	5,5 x 8,7	5,5 x 8,7	

Set EY-3

20x EY-3m | 20x EY-3f



Set EY-1

20x EY-1m | 20x EY-1f



3.000 - 8.000 min⁻¹
max. 20.000 min⁻¹



40 pieces / kosov

Item / Artikel:	Set EY-3	EY-3m	EY-3f	Set EY-1	EY-1m	EY-1f	RA
REF:	9031	3273	3373	9030	3271	3371	
Dimensions / Dimenzija (mm):	—	6,2 x 6,4	6,2 x 6,4	—	4 x 6,9	4 x 6,9	

ECOCOMP

Composite · Kompozit




POLISHING SYSTEM | POLIRNI SISTEM

- For efficient polishing of composite surfaces
- Silicon carbide grit
- Hard bonding for high durability
- Za efektivno poliranje kompozitnih površin
- Silicijeva karbidna zrnca
- Trdna vezava za visoko vzdržljivost

Set RA 210 REF 9071

EC11m | EC7m | EC3Km | EC18m
EC11f | EC7f | EC3Kf | EC18f



1

MEDIUM Pre-polishing
SREDNJI Predpoliranje

2

FINE Polishing
FINI Poliranje

2

PROCESSING STEPS
STOPENJSKI PROCES



7.000 - 10.000 min⁻¹
max. 20.000 min⁻¹



10, 100 pieces / kosov

RA	EC11m	EC3Km	EC7m	EC11f	EC3Kf	EC7f	Item / Artikel:
	5181	5173	5177	5281	5273	5277	REF:
	3 x 7,5	6,5 x 9,5	5 x 10	3 x 7,5	6,5 x 9,5	5 x 10	Dimensions / Dimenzija (mm):

COMPOSOFT

Composite · Kompozit

POLISHING SYSTEM | POLIRNI SISTEM

- For smoothing, pre-polishing and polishing composite
- Aluminium oxide grit
- Soft silicone bonding adapts to the composite surface
- Za glajenje, predpoliranje in poliranje kompozita
- Alumijevo oksidna zrnca
- Meha silikonska vezava se prilagaja površini kompozita

Set RA 129 REF 9089

CS11g | CS2Bg | CS7g
CS11m | CS2Bm | CS7m
CS11f | CS2Bf | CS7f



1

COARSE Smoothing
GROBI Glajenje

2

MEDIUM Pre-polishing
SREDNJI Predpoliranje

3

FINE Polishing
FINI Poliranje

3

PROCESSING STEPS
STOPENJSKI PROCES



5.000 - 8.000 min⁻¹
max. 20.000 min⁻¹



10, 100 pieces / kosov

RA	CS11g	CS2Bg	CS7g	CS11m	CS2Bm	CS7m	CS11f	CS2Bf	CS7f	Item / Artikel:
	2581	2592	2577	2681	2692	2677	2781	2792	2777	REF:
	3 x 7,5	5 x 10	7 x 10	3 x 7,5	5 x 10	7 x 10	3 x 7,5	5 x 10	7 x 10	Dimensions / Dimenzija (mm):

DIAPOL®

Feldspar Ceramic, Lithium Disilicate, Zirconia
Leucitna Keramika, Litijev Disilikat, Cirkonija



DIAMOND IMPREGNATED POLISHING SYSTEM
Z DIAMANTI IMPREGNIRAN POLIRNI SISTEM

- 3-step polishing system for all porcelain materials
- Diamond grit for highest efficiency
- Surface corrections without additional glazing
- 3 stopenjski polirni sistem za vse porcelanske materiale
- Diamantni delci za visoko učinkovitost
- Površinski popravki brez dodatne glazure

Set RA 308 REF 9060

OF-Dg | OF-Dmf | OF-D

Set RA 306 REF 9088

DT-W14Dg | DT-W14Dmf | DT-W14D

Set RA 307 REF 9061

W16Dg | OF-Dg | W17Dg | W18Dg
W16Dmf | OF-Dmf | W17Dmf | W18Dmf
W16D | OF-D | W17D | W18D

3 PROCESSING STEPS STOPENJSKI PROCES

1 COARSE Smoothing
GROBI Glajenje

2 MEDIUM Pre-polishing
SREDNJI Predpoliranje

3 FINE High-shine polishing
FINI Poliranje visoki sijaj

7.000 - 12.000 min⁻¹
* 3.000 - 8.000 min⁻¹
max. 20.000 min⁻¹

1, 10 piece / kosov



Item / Artikel:	DT-W14Dg	OF-Dg	W11Dg	W16Dg	W17Dg	RA
REF:	7084	7075	7081	7086	7087	
Dimensions / Dimenzija (mm):	14 x 1,6	5,5 x 8,8	3 x 7,5	4 x 10	6 x 7,5	

7.000 - 12.000 min⁻¹
* 3.000 - 8.000 min⁻¹
max. 20.000 min⁻¹

1, 10 piece / kosov



Item / Artikel:	DT-W14Dmf	OF-Dmf	W11Dmf	W16Dmf	W17Dmf	RA
REF:	7184	7175	7181	7186	7187	
Dimensions / Dimenzija (mm):	14 x 1,6	5,5 x 8,8	3 x 7,5	4 x 10	6 x 7,5	

7.000 - 12.000 min⁻¹
* 3.000 - 8.000 min⁻¹
max. 20.000 min⁻¹

1, 10 piece / kosov



Item / Artikel:	DT-W14D	OF-D	W11D	W16D	W17D	RA
REF:	7284	7275	7281	7286	7287	
Dimensions / Dimenzija (mm):	14 x 1,6	5,5 x 8,8	3 x 7,5	4 x 10	6 x 7,5	

DIACERA

Zirconia · Cirkonija



DIAMOND IMPREGNATED POLISHING SYSTEM
Z DIAMANTI IMPREGNIRAN POLIRNI SISTEM

- Especially suitable for full contour zirconia restorations
- High diamond concentration for more efficiency
- Avoids repeated glazing after retouching
- Posebej primerno za polne cirkonske restavracije
- Visoka vsebnost diamantov za večjo učinkovitost
- Ni potrebe za dodatno glazuro po popravkih

NEW

Set RA 326 REF 9062

OF-DCmf | OF-DC



Set RA 323 REF 9091

DT-W14DCmf | DT-W14DC

NEW

Set RA 325 REF 9063

W16DCmf | OF-DCmf | W17DCmf | W18DCmf
W16DC | OF-DC | W17DC | W18DC

+ DIASYN[®] plus (Grinding · Brusilni)
DYP-W13g



MEDIUM Pre-polishing
SREDNJI Predpoliranje



FINE High-shine polishing
FINI Poliranje visoki sijaj



2 PROCESSING STEPS
STOPENJSKI PROCES



7.000 - 10.000 min⁻¹
* 3.000 - 8.000 min⁻¹
max. 20.000 min⁻¹
1, 10 piece / kosov

RA	DT-W14DCmf	OF-DCmf	W11DCmf	W16DCmf	W17DCmf	Item / Artikel:
	7684	7675	7681	7686	7687	REF:
	14 x 1,6	5,5 x 8,8	3 x 7,5	4 x 10	6 x 7,5	Dimensions / Dimenzija (mm):



7.000 - 10.000 min⁻¹
* 3.000 - 8.000 min⁻¹
max. 20.000 min⁻¹
1, 10 piece / kosov

RA	DT-W14DC	OF-DC	W11DC	W16DC	W17DC	Item / Artikel:
	7784	7775	7781	7786	7787	REF:
	14 x 1,6	5,5 x 8,8	3 x 7,5	4 x 10	6 x 7,5	Dimensions / Dimenzija (mm):

DIAPRO

Lithium Disilicate · Litijev Disilkat



DIAMOND IMPREGNATED POLISHING SYSTEM
Z DIAMANTI IMPREGNIRAN POLIRNI SISTEM

- Especially suitable for silicate porcelain
- Diamond grit for highest efficiency
- Polishing result matches glazing
- Posebej primerno za silikatni porcelan
- Diamantni delci za najvišjo učinkovitost
- Polni rezultati kot pri glaziranju

NEW

Set RA 366 REF 9064

OF-DPmf | OF-DP



Set RA 363 REF 9092

DT-W14DPmf | DT-W14DP



NEW

Set RA 365 REF 9065

W16DPmf | OF-DPmf | W17DPmf | W18DPmf
W16DP | OF-DP | W17DP | W18DP

+ DIASYN[®] plus (Grinding · Brusilni)
DYP-W13m



2

PROCESSING STEPS
STOPENJSKI PROCES



MEDIUM Pre-polishing
SREDNJI Predpoliranje



FINE High-shine polishing
FINI Poliranje visoki sijaj



7.000 - 10.000 min⁻¹
* 3.000 - 8.000 min⁻¹
max. 20.000 min⁻¹



1, 10 piece / kosov



Item / Artikel:

DT-W14DPmf

OF-DPmf

W11DPmf

W16DPmf

W17DPmf

RA

REF:

7884

7875

7881

7886

7887

Dimensions / Dimenzija (mm):

14 x 1,6

5,5 x 8,8

3 x 7,5

4 x 10

6 x 7,5

NEW



7.000 - 10.000 min⁻¹
* 3.000 - 8.000 min⁻¹
max. 20.000 min⁻¹



1, 10 piece / kosov

Item / Artikel:

DT-W14DP

OF-DP

W11DP

W16DP

W17DP

RA

REF:

7984

7975

7981

7986

7987

Dimensions / Dimenzija (mm):

14 x 1,6

5,5 x 8,8

3 x 7,5

4 x 10

6 x 7,5

POLISHING SYSTEM | POLIRNI SISTEM

- For finishing and polishing porcelain restorations
- Silicone bonding adapts to the surface being processed
- Silicon carbide abrasive
- Za končno obdelavo in poliranje porcelanskih restavracij
- Silikonska vezava se prilagaja obdelovani površini
- Silicijevi karbidni abrazivi

Set RA 160 REF 9076

W2BVK | W7VK
W2BNK | W7NK
W2BSK | W7SK



1 COARSE Smoothing GROBI Glajenje



2 MEDIUM Pre-polishing SREDNJI Predpoliranje



3 FINE Polishing FINI Poliranje



3 PROCESSING STEPS STOPENJSKI PROCES

7.000 - 15.000 min⁻¹
max. 25.000 min⁻¹

10, 100 pieces / kosov

RA	W2BVK	W7VK	W2BNK	W7NK	W2BSK	W7SK	Item / Artikel:
	1592	1577	1692	1677	1792	1777	REF:
	5 x 10	7 x 10	5 x 10	7 x 10	5 x 10	7 x 10	Dimensions / Dimenzija (mm):

UNIVERSAL

Composite, Feldspar Ceramic, Amalgam, Gold
Kompozit, Leucitna Keramika, Amalgam, Zlato

POLISHING SYSTEM | POLIRNI SISTEM

- For smoothing all dental materials
- Soft silicone bonding adapts to the surface being processed
- Silicon carbide abrasive
- Za glajenje vseh zobozdravstvenih materialov
- Mehka vezava silikona se prilagaja obdelovani površini
- Silicijevi karbidni abrazivi

Set RA 101 REF 9070

W11 | W1 | W2 | W2B | W5 | W7 | W8 | W9



1 MEDIUM Smoothing SREDNJI Predpoliranje



3, 10, 100 pieces / kosov

1 PROCESSING STEP STOPENJSKI PROCES

7.000 - 10.000 min⁻¹
max. 20.000 min⁻¹

10, 100 pieces / kosov

RA	K11	76RA	W11	W1	W5	W7	Item / Artikel:
	1011	9576	1081	1071	1075	1077	REF:
	11 x 6	5	3 x 7,5	5,5 x 15	9 x 8,5	7 x 10	Dimensions / Dimenzija (mm):

EVEFLEX

Amalgam, Gold · Amalgam, Zlato

POLISHING SYSTEM | POLIRNI SISTEM

- For polishing amalgam and gold
- 3-step polishing system
- Silicon carbide grit
- Za poliranje amalgama in zlata
- 3 stopenjski polirni sistem
- Silicijevi karbidni abrazivi

Set RA 240 REF 9090

ET-AF14 | ET-BV14 | ET-GH14



Set RA 209 REF 9075

AF11 | AF7 | AF3K
BV11 | BV7 | BV3K
GH11 | GH7 | GH3K



3 PROCESSING STEPS STOPENJSKI PROCES

1 COARSE Smoothing
GROBI Glajenje

2 MEDIUM Pre-polishing
SREDNJI Predpoliranje

3 FINE Polishing
FINI Poliranje

7.000 - 12.000 min⁻¹
max. 20.000 min⁻¹

10, 100 pieces / kosov



Item / Artikel:	ET-AF14	AF11	AF3K	AF7	RA
REF:	4684	4681	4673	4677	
Dimensions / Dimenzija (mm):	14 x 1,6	3 x 7,5	6,5 x 9,5	5 x 10	

7.000 - 12.000 min⁻¹
max. 20.000 min⁻¹

10, 100 pieces / kosov



Item / Artikel:	ET-BV14	BV11	BV3K	BV7	RA
REF:	4784	4781	4773	4777	
Dimensions / Dimenzija (mm):	14 x 1,6	3 x 7,5	6,5 x 9,5	5 x 10	

7.000 - 12.000 min⁻¹
max. 20.000 min⁻¹

10, 100 pieces / kosov



Item / Artikel:	ET-GH14	GH11	GH3K	GH7	RA
REF:	4884	4881	4873	4877	
Dimensions / Dimenzija (mm):	14 x 1,6	3 x 7,5	6,5 x 9,5	5 x 10	

FLEXI-DIA



Feldspar Ceramic, Lithium Disilicate, Zirconia, Hybrid Ceramic
Leucitna Keramika, Litijev Disilikat, Cirkonija, Hibridna Keramika

DIAMOND COATED POLISHING SYSTEM
Z DIAMANTI IMPREGNIRAN POLIRNI SISTEM

- Flexible polishing discs with diamond coating (single side)
- For use in a pulling or pushing motion
- Extra-thin plastic foil
- Prilagodljivi polirni diski z diamantnim nanosom (enostranski)
- Uporabni pri vlečnih ali porivnih gibih
- Izjemno tanek kolot



3, 10, 100 pieces / kosov



7.000 - 10.000 min⁻¹
max. 20.000 min⁻¹

100 pieces / kosov
available with mandrel
na voljo z mandrelo

RA	FD-14Dxg	FD-14Dg	FD-14Dmf	FD-14D	95RA	Item / Artikel:
	9146	9147	9148	9149	9595	REF:
	14 x 0,21	14 x 0,17	14 x 0,14	14 x 0,11	-	Dimensions / Dimenzija (mm):

FLEXI-D

Composite, Amalgam, Gold
Kompozit, Amalgam, Zlato

POLISHING SYSTEM | POLIRNI SISTEM

- Flexible polishing discs coated with aluminium oxide (single side)
- For surface preparation and interproximal polishing
- For use in a pulling or pushing motion
- Fleksibilni polirni diski prevlečeni z aluminijevim oksidom (enostranski)
- Za površinsko preparacijo in interproksimalno poliranje
- Uporabni z vlečnimi ali porivnimi gibi



3, 10, 100 pieces / kosov



7.000 - 10.000 min⁻¹
max. 20.000 min⁻¹

100 pieces / kosov
available with mandrel
na voljo z mandrelo

RA	FD-14g	FD-14m	FD-14f	FD-14sf	95RA	Item / Artikel:
	9141	9142	9143	9144	9595	REF:
	14 x 0,21	14 x 0,17	14 x 0,14	14 x 0,11	-	Dimensions / Dimenzija (mm):

EVEBRUSH

Composite · Kompozit

POLISHING SYSTEM | POLIRNI SISTEM

- For high-shine polishing on composite
- For polishing occlusal fillings
- Fibre with silicon carbide grit
- Application without polishing paste
- Heat resistant and autoclavable polyamide fibre

- Za poliranje kompozitov do visokega sijaja
- Za poliranje okluzijskih polnitev
- Vlakna z delci silicijevega karbida
- Uporaba brez polirne paste
- Avtoklavirna poliamidna vlakna odporna na vročino



FINE: High-shine polishing
HARDNESS: hard

FINI: Poliranje visoki sijaj
TRDOTA: trdi



max. 10.000 min⁻¹



1, 10 piece / kosov



Item / Artikel:	WT1	WZ1	WS1	RA
REF:	9491	9492	9493	
Dimensions / Dimenzija (mm):	6 x 5,5	4 x 6,5	4 x 6,5	

BRACKET *plus*

Debonding · Odstranjevanje

POLISHING SYSTEM | POLIRNI SISTEM

- For removal of remnants after debonding
- Does not harm the enamel
- PROPHYCLEAN is recommended for removal of hard bracket remnants

- Odstranjevanje ostankov po odstranitvi orto nosilcev
- Ne poškoduje sklenine
- PROPHYCLEAN se priporoča za odstranjevanje trdih ostankov nosilcev

Set RA 431 REF 9097

BRP11 | BRP7 | BRP3 | BRP6

+ PROPHYCLEAN
(Selective removal · Selektivno odstranjevanje)
PC-100



5.000 - 7.000 min⁻¹
max. 20.000 min⁻¹



10, 100 pieces / kosov



Item / Artikel:	BRP11	BRP7	BRP3	BRP6	RA
REF:	5681	5677	5673	5676	
Dimensions / Dimenzija (mm):	3 x 7,5	5 x 10	5,5 x 8,5	10 x 0,7	


PROPHYCLEAN



Tooth Cleaning, Debonding
Čiščenje Zoba, Odstranjevanje

POLISHING SYSTEM | POLIRNI SISTEM

- Removes calculus and stains
- Special abrasive prevents damage to the enamel
- For mechanical prophylactic treatment
- Fibre free
- Blue option for heavy stains
- Odstranjuje zobni kamen in madeže
- Poseben abraziv preprečuje poškodbe sklenine
- Za mehansko profilaktično obdelavo
- Brez vlaken
- Modra izbor za močno obarvanje



- 5.000 - 10.000 min⁻¹
max. 20.000 min⁻¹
- 1, 10 piece / kosov

RA	PC-100g	PC-102g	PC-103g	PC-100	PC-102	PC-103	Item / Artikel:
	8700	8702	8703	8800	8802	8803	REF:
	2,35 x 9	2,35 x 7	2,35 x 7,5	2,35 x 9	2,35 x 7	2,35 x 7,5	Dimensions / Dimenzija (mm):

DIASTRIPS

Composite, Feldspar Ceramic, Lithium Disilicate, Zirconia, Hybrid Ceramic, NEM Alloys, Titanium
Kompozit, Leucitna Keramika, Litijev Disilikat, Cirkonija, Hibridna Keramika, NEM Kovine, Titan

DIAMOND COATED SEPARATING STRIPS
Z DIAMANTI IMPREGNIRAN POLIRNI SISTEM

- For manual surface processing
- Interdental application
- Perforation removes the cut material
- Za ročno obdelavo površine
- Interdentalna uporaba
- Perforacija odstranjuje višek materiala

Set PDL4 REF 9014
PDL-45/4 | PDL-35/4 | PDL-10/4

Set PDL2,5 REF 9012
PDL-45/2,5 | PDL-35/2,5 | PDL-10/2,5



	PDL-45/4	PDL-45/2,5	PDL-35/4	PDL-35/2,5	PDL-10/4	PDL-10/2,5	Item / Artikel:
	9101	9106	9102	9107	9103	9108	REF:
	4 x 150	2,5 x 150	4 x 150	2,5 x 150	4 x 150	2,5 x 150	Dimensions / Dimenzija (mm):
	0,11	0,11	0,10	0,10	0,07	0,07	Thickness / Debelina (mm):

- 9 pieces / kosov



Download brochure
Prejmite bršuro

MADE IN GERMANY

Perfect Surfaces.

EVE Ernst Vetter GmbH
Neureutstr. 6
75210 Kelttern, Germany

+49 (0) 72 31 97 77 -0
info@eve-rotary.com

www.eve-rotary.com


SCAN ME!

**Interested in
EVE products?**

Scan the QR code and find an
EVE distributor in your country.

**Vas zanimajo
EVE izdelki?**

Skenirajte kodo QR in poiščite
distributerja EVE v svoji državi.



The specifications and appearance of the products are subject to change without notice. Printed colour can be slightly different from actual colour. Changes, variations and errors reserved. Notes: www.eve-rotary.com/downloads

Specifikacije in videz izdelkov se lahko spremenijo brez obvestila. Natisnjena barva se lahko nekoliko razlikuje od dejanske barve. Pridržujemo si pravico sprememb, variacij in napak. Opombe: www.eve-rotary.com/downloads

B005_Ed11.1-A_EN_SL