

# SPECIAL OFFER | Clinic

From 1 February to 31 May, 2025



## icanCLAVE 3F Autoclave



icanCLAVE 3F is a super-fast autoclave, especially for the sterilisation of unwrapped instruments, which greatly improves the rotation rate of surgical and the convenience of small operations such as dentistry.

- ✓ Printer port.
- ✓ USB port.
- ✓ Internal memory.
- ✓ Barcode label printer port.
- ✓ Carrying case.




Official technical service  
**icanCLAVE**  
in Spain

REF. 100367

LAUNCH OFFER  
**1.995 €**


### Specifications

Voltage, Frequency: 220/110V, 50/60Hz • Chamber volume: Aprox. 3L  
Power: 2300 W • Number of trays: 1 • Tray material: Aluminium  
Tray dimensions: 178 x 285 mm  
Chamber dimension: 192(W) x 45(H) x 285(D)



### REGULATIONS

- ✓ 97/23/EEC Pressure Equipment Directive.
  - ✓ Regulation (EU) 2017/745 on medical devices (MDR)
  - ✓ EN13060 Small steam sterilizers.
  - ✓ EN61326, EN 61000 Electromagnetic compatibility.
  - ✓ EN61010-1, UL 61010-1 Safety requirements.
  - ✓ EN61010-2-040, IEC 61010-2-040 Specific requirements for steam sterilizers used for treatment of medical materials.
  - ✓ EN 13485:2012/AC:2012
  - ✓ ISO 13485:2003
- Official technical service ICANCLAVE in Spain.



## 18-D Desktop autoclave class B (18 liters)



- ✓ High quality internal components.
- ✓ ODE valves, ULUAK vacuum pump, ULKA water pump.
- ✓ Ultra quiet (70 dB).
- ✓ Triple fractionated vacuum system with an independent steam generator, which guarantees fast and safe sterilization cycles.
- ✓ High precision control system for perfect sterilization results (EN13060).
- ✓ Auto-diagnosis with error codes.
- ✓ Fully automatic with pre-programmed sterilization cycles.
- ✓ Chamber and door in stainless steel, with electro-polished finish.
- ✓ Open type water tank, which allows easy filling and cleaning.
- ✓ With printer.



Official technical service  
**icanCLAVE**  
in Spain

REF. 100365

**2.852 €**

### REGULATIONS

- ✓ 97/23/EEC Pressure Equipment Directive.
- ✓ Regulation (EU) 2017/745 on medical devices (MDR)
- ✓ EN13060 Small steam sterilizers.
- ✓ EN61326, EN 61000 Electromagnetic compatibility.
- ✓ EN61010-1, UL 61010-1 Safety requirements.
- ✓ EN61010-2-040, IEC 61010-2-040 Specific requirements for steam sterilizers used for treatment of medical materials.
- ✓ EN 13485:2012/AC:2012
- ✓ ISO 13485:2003



### Specifications

Chamber dimensions 247 x 350 mm  
Nominal power 1750 V  
Operation temperature 5-40°C  
Operation relative humidity Máx. 80 %  
Operation atmospheric pressure 76-106 kPa  
Noise < 70 dB  
Circuit breaker F16A/400 V  
Sterilization temperatures 121/134°C  
Distilled water tank capacity 2.5 l

## Suction unit Sirocco S



Incorporates two stages of filtration: a first circular HEPA (High Efficiency Particulate Array) filter and a second internal HEPA filter. The unit also features an external micromotor detection system and seven different power levels.

Together with the finishing box, they form an extremely useful and reliable partnership for a clinic environment, allowing for minor polishing and reworking operations.

- ✓ Powerful induction motor.
- ✓ Intuitive touch keypad.
- ✓ Easy to use and change filter.
- ✓ Large screen and access.
- ✓ Exceptionally quiet.
- ✓ LED lamp for perfect illumination.
- ✓ Great working ergonomics.



Sirocco S + Finishing box

REF. 080534AC

**699 €**

### Specifications

Box dimensions 365 x 290 x 290 mm • Weight 5.5 kg  
Power 350 W • Noise 55 dB

## Biological espore incubators

Simple and easy to use, they are essential to carry out the necessary tests that confirm the proper performance of the sterilization devices.

Both have four cameras to be able to perform as many verifications at the same time. The Elite model's electronic control also alerts

of anomalies such as power outages to avoid incorrect results.

Essentials to confirm the proper performance of the autoclave.



Elite model  
Adjustable temperature  
and digital display

REF. 100370

**300 €**



Practic model  
Functional  
and economical

REF. 100375

**129 €**

### Specifications

Power 25 W  
Temperature range 57.5 ± 2.5°C / 35 ± 2.5°C

### Specifications

Power 10 W  
Temperature range 57.5°C

# SPECIAL OFFER | Clinic

From 1 February to 31 May, 2025



MESTRA

## Pressure former PlastPress

The most compact positive pressure machine on the market.

Greater accuracy in the forming process and lower noise compared to vacuum machines.

Practical and versatile, it accepts round and square sheets up to 6 mm thickness.

It must be connected directly to a compressor (between 3 and 8 bar).

Sheets  $\square$  from 125 to 131 mm  
Sheets  $\circ$  from 125 to 134 mm



NEW

REF. 080544

941 €

### Specifications

Height 270 mm  
Width 235 mm  
Depth 265 mm  
Weight 7.5 kg  
Power 450 W

## Vacuum forming machine

Thermoforming machine for square and round sheets.

Sheets  $\square$  from 125 to 134 mm  
Sheets  $\circ$  from 125 to 134 mm  
(adapter available for  $\circ$  120, not included)



REF. 080542

651 €

### Specifications

Height 310 mm  
Width 275 mm  
Depth 260 mm  
Weight 4.6 kg  
Power 925W

NEW

## Thermoforming sheets

Double packaging (inside-outside).  
For all types of use.

CE Medical device



REFERENCE	HARDNESS	APPEARANCE	Dim. 126 x 126 mm thickness	QTY.	PRICE
080563	Hard	Transparent	0.5 mm	25 u.	11 €
080564	Hard	Transparent	0.75 mm	25 u.	14,56 €
080565	Hard	Transparent	1.0 mm	25 u.	14,56 €
080566	Hard	Transparent	1.5 mm	25 u.	18,30 €
080567	Hard	Transparent	2.0 mm	25 u.	27 €
080568	Soft	Translucent	0.5 mm	25 u.	25 €
080569	Soft	Translucent	1.0 mm	12 u.	6,98 €
080570	Soft	Translucent	1.5 mm	12 u.	9,27 €
080571	Soft	Translucent	2.0 mm	12 u.	13,97 €
080572	Soft	Translucent	3.0 mm	12 u.	17,45 €

## Thermoforming sheets

CE Medical device



Standard square and round thermoforming sheets. Made of high quality plastic with CE marking certification as medical device.

There is no soot, nor fumes during the heating. As they do not contain recycled materials, they do not tear while subjected to vacuum.

Available in a wide variety of thicknesses and types. Supplied in metalized package and vacuum closed.

From 4 €

REF.	TYPE	HARDNESS	APPEARANCE	Dim. 126x126 mm thickness	QTY.	OFFER PRICE
080574	Crystal	Hard	Transparent	0.3 mm	10	4€
080575	Crystal	Hard	Transparent	0.75 mm	5	4€
080576	Crystal	Hard	Transparent	1.0 mm	5	4€
080577	Crystal	Hard	Transparent	1.5 mm	5	6€
080578	Crystal	Hard	Transparent	2.0 mm	5	8€
080579	Milky	Hard	Translucent	0.6 mm	10	4€
REF.	TYPE	HARDNESS	APPEARANCE	Dim. 131x131 mm thickness	QTY.	OFFER PRICE
080580	Soft	Soft	Translucent	1.0 mm	10	7€
080581	Soft	Soft	Translucent	1.5 mm	10	10€
080582	Soft	Soft	Translucent	2.0 mm	10	12€
080583	Soft	Soft	Translucent	3.0 mm	5	9€
080584	Soft	Soft	Translucent	4.0 mm	5	13€
REF.	TYPE	HARDNESS	APPEARANCE	Dim. $\phi$ 125 mm thickness	QTY.	OFFER PRICE
080574R	Crystal	Hard	Transparent	0.3 mm	10	4€
080575R	Crystal	Hard	Transparent	0.75 mm	5	4€
080576R	Crystal	Hard	Transparent	1.0 mm	5	5€
080577R	Crystal	Hard	Transparent	1.5 mm	5	7€
080579R	Milky	Hard	Translucent	0.6 mm	10	5€
080580R	Soft	Soft	Translucent	1.0 mm	10	5€
080581R	Soft	Soft	Translucent	1.5 mm	10	9€
080582R	Soft	Soft	Translucent	2.0 mm	10	9€
080583R	Soft	Soft	Translucent	3.0 mm	5	6€
080584R	Soft	Soft	Translucent	4.0 mm	5	11€

## Pressure former



The strong pressurization allows precise molding, regardless of the cut of the model, as well as various thicknesses (0.5 to 4.0 mm). Square and round sheets (127, 125, 120 mm) can be used.

### Specifications

Height 545 mm Width 410 mm  
Depth 260 mm Weight 15 kg  
Power 450 W Pressure 4-6 bar

REF. 080543

2.422 €

## Drying oven for thermoforming sheets



A fan circulates the hot air producing a uniform drying:

- ✓ Draws out moisture and keeps dry thermoforming sheets, to avoid the formation of bubbles.
- ✓ Drying of instruments plaster models, etc.
- ✓ 6 stainless steel trays.
- ✓ Stainless steel housing.
- ✓ 24 hour timer.

REF. 080118

268 €

### Specifications

Height 310 mm Width 340 mm Depth 450 mm  
Weight 8.3 kg Power 550-650 W Temperature 35-75 °C



# SPECIAL OFFER | Clinic

From 1 February to 31 May, 2025



MESTRA

## Automatic bag sealer



NEW

REF. 080435

2.170 €

200 mm  
240 mm  
635 mm

Optional roller  
REF. 080436

17 kg

## Bag sealer



Bags not included

### Specifications

Height 210 mm  
Width 430 mm  
Depth 335 mm  
Weight 6.5 kg  
Power 100 W

REF. 080433

SUPER OFFER  
265 €

- ✓ Sealing length 305 mm
- ✓ Sealing width 12 mm

## Bags for sealer self-sealing

New bags 5x9 cm  
The smallest on the market!



Class I medical device  
(MDR 2017/745)  
Conforming to standards  
ISO 11607-1:2017 and  
EN 868-5:2009

NEW


### Specifications

Ref. 080448: 5 x 24 cm	Ref. 080449: 7 x 26 cm
Ref. 080450: 9 x 26 cm	Ref. 080451: 5 x 9 cm
Ref. 080452: 14 x 28 cm	Ref. 080454: 19 x 36 cm
Ref. 080456: 25 x 40 cm	Ref. 080457: 30 x 40 cm
Ref. 080458: 30 x 45 cm	Ref. 080459: 5.7 x 13 cm

200 units  
per package

CHECK  
PRICES

## Coils for sealers



### Specifications

Ref. 080441: 5 cm x 200 m  
Ref. 080442: 7.5 cm x 200 m  
Ref. 080443: 10 cm x 200 m  
Ref. 080444: 15 cm x 200 m  
Ref. 080445: 20 cm x 200 m  
Ref. 080446: 25 cm x 200 m  
Ref. 080447: 30 cm x 200 m

Class I medical device  
(MDR 2017/745)  
Conforming to standards  
ISO 11607-1:2017 and EN 868-5:2009

CHECK  
PRICES

## Barrier film

This adhesive film offers a practical and effective solution for dental environments, providing a layer of protection for instruments and materials.

Each roll contains 1200 pre-cut 10 x 15 cm adhesives.  
Blue colour.



NEW

REF. 080439

LAUNCH OFFER  
10 €

### Specifications

1200 pre-cut adhesives.  
Dimensions 10 x 15 cm

## Germicide box for burs



Stainless steel

REF. 070530

42 €

REF. 070531

8 €

### Specifications

Height 80 mm  
Diameter 90 mm  
Weight 0.12 kg

### Specifications

Height 35 mm  
Diameter 70 mm  
Weight 0.12 kg

## Tank for distilled water

Made of stainless steel, it features a pressure regulator and a double safety valve

Inox

MADE IN  
MESTRA  
SPAIN  
DESIGN



### Specifications

Height 680 mm  
Width (with pressure regulator) 255 mm  
Diameter 225 mm  
Weight 4.7 kg  
Capacity 18 l

REF. 080620

307 €

## Distiller decalcifier

### Specifications

Height 385 mm  
Width 250 mm  
Depth 240 mm  
Weight 3.4 kg  
Capacity 4 l

REF. 100210

230 €



# SPECIAL OFFER | Clinic

From 1 February to 31 May, 2025



MESTRA

## Ultrasonics bath 2.8, 3.4 and 5.8 liters

- ✓ Temperature.
- ✓ Drain.
- ✓ Metal basket.
- ✓ Made of Stainless Steel.

Inox



5.8 L  
REF. 100152

452 €

Inner dimensions:  
Height 155 mm  
Width 298 mm  
Depth 152 mm



3.4 L  
REF. 100151

421 €

Inner dimensions:  
Height 100 mm  
Width 298 mm  
Depth 152 mm

NEW



2.8 L  
REF. 100150

328 €

Inner dimensions:  
Height 100 mm  
Width 238 mm  
Depth 140 mm

## Calypso ultrasonic bath 4 liters



- ✓ With heating
- ✓ With drain valve
- ✓ Includes crystal glass
- ✓ Includes cutter cleaning pliers (ref. 070520)
- ✓ Includes stainless steel basket

REF. 080422

1.124 €

### Specifications

Height 220 mm • Width 425 mm • Depth 270 mm  
Weight 4.5 kg • Capacity 4 l • Frequency 40 kHz

## Ultrasonic bath 180 ml with UV-light

NEW




REF. 100162

LAUNCH OFFER  
47 €

### Specifications

Height 122 mm • Width 122 mm • Depth 104 mm  
Weight 2,8 kg • Capacity 180 ml • Frequency 42 kHz  
Cuve Ø 87\*38 mm • Timer 5/10 min  
Tension 220 50/60 Hz • Power 28W

## Ultrasonic bath For denture cleaning



75 mm  
85 mm  
90 mm  
0.18 kg

### Specifications

Power supply 2 batteries AA

REF. 100160

\*The battery is not included.

CHECK  
PRICES

## Ultrasonic bath 3, 2.5 and 1.4 liters

NEW



REF. 100140

262 €



REF. 100142

170 €



REF. 100141

170 €

REF.	Height	Width	Depth	Weight	Capacity
100140	204 mm	314mm	243 mm	3,1 kg	2,5 L
100141	162 mm	225 mm	240 mm	1,9 kg	1,4 L
100142	234 mm	385 mm	193 mm	3,6 kg	3 L



# SPECIAL OFFER | Clinic

From 1 February to 31 May, 2025



MESTRA

## Alginate mixer

- ✓ Transportable.
- ✓ Ergonomic.
- ✓ Two speeds.
- ✓ Spatule included.

CE Medical device  
MDR 2017/745



Gift - One cup  
valued in 27,69 €  
Ref. 080515-02



### Specifications

Height 230 mm  
Width 175 mm  
Depth 250 mm  
Weight 3 kg

REF. 080515

403 €

## Alginate automatic mixer

- ✓ Comfortable automatic mixing.
- ✓ Quick mixing obtention in a few seconds.
- ✓ Easy cleaning of all components.
- ✓ Max mix 80 grs.



Gift - Sealed alginate  
container valued  
at 93,20 €  
Ref. 100189



305 mm  
265 mm  
205 mm  
17 kg

REF. 100191

819 €

## Apex locator

- ✓ Automatic calibration.
- ✓ LCD display clear and accurate image.
- ✓ Different colours.
- ✓ Sterilizable accessories.
- ✓ Rechargeable battery.



REF. 100330

SUPER OFFER  
291 €

110 mm  
80 mm  
138 mm  
0.4 kg

## Intraoral camera (autofocus)

- ✓ High resolution with 2000 000 pixels.
- ✓ Transfer of images with WIFI.
- ✓ 15 inches LCD screen.
- ✓ Autofocus function.



REF. 100335

749 €

360 mm  
150 mm  
308 mm  
3 kg

## LED C curing light

- ✓ It has a desktop base/charger.
- ✓ Without wires.

CE  
0197

### Specifications

Dimensions 263 x Ø 23 mm  
Weight 120 gr  
Power supply 100-240 V 50/60 Hz  
Light output 850-1000 mW/cm²



REF. 100320

SUPER OFFER  
190 €

## Motor for endodontics

CE  
0120

### Specifications

Height 215 mm  
Width 85 mm  
Depth 120 mm  
Weight 0.3 kg  
Speed 100-600 rpm  
Torque 1.4 N. cm  
Power supply Batteries



REF. 100350

SUPER OFFER  
612 €

## Tabletop ultrasonic scaler

CE  
0197

80 mm  
134 mm  
200 mm  
0.66 kg



D3

REF. 100313

250 €

2 Functions  
Scaling - Periodontics  
Satelec type handpiece  
Detachable and autoclavable with 5 tips.

81 mm  
159 mm  
188 mm  
0.65 kg



D5

REF. 100316

SUPER OFFER  
300 €

3 Functions  
Scaling - Periodontics - Endodontics  
Satelec type handpiece  
Detachable and autoclavable with 6 tips.

230 mm  
190 mm  
240 mm  
1.23 kg



D7

REF. 100318

SUPER OFFER  
622 €

3 Functions  
Scaling - Periodontics - Endodontics  
Satelec type handpiece. Detachable and autoclavable at 135° with 8 tips. With liquid reservoir or drug treatment.

# SPECIAL OFFER | Clinic


From 1 February to 31 May, 2025



MESTRA

## Whitening Lamp

- ✓ Wide spectrum of illumination (from 420 to 472 nm)
- ✓ High power (up to 2000 mW/cm<sup>2</sup>)
- ✓ Easily adjustable in height and orientation.
- ✓ Its digital display shows the power, times, etc. of the whitening process.
- ✓ Includes remote control 20 disposable mouth splints, two protective glasses and a lip retractor.



### Specifications

Height 1120 mm
Width 820 mm
Depth 280 mm
Weight 10 kg

REF. 100340

**SUPER OFFER**  
**490 €**

## Whitening Lamp

- ✓ Wide spectrum of illumination.
- ✓ 3 different types of LED light sources: blue, red, (red + blue).
- ✓ 5 blue lights from 430 to 490 nm.
- ✓ 2 high power red lights 610 to 650 nm.
- ✓ Light output power up to 2327 mW/cm<sup>2</sup>
- ✓ 8 purple lights with sterilization effect.
- ✓ Made of metallic material, in appearance elegant and easy to assemble.
- ✓ Adjustable bleaching time.



### Specifications

Height 1700 mm
Width 560 mm
Depth 560 mm
Weight 19 kg

REF. 100341

**691 €**

## Mini sandblaster for extraoral use

- ✓ Particle size: 50 to 100 microns.
- ✓ Essential for brackets cleaning.
- ✓ Adequacy of micro retentions and crowns.
- ✓ Adjustment of interproximal contacts.
- ✓ Occlusal sandBlasting.
- ✓ Inlay-Onlay.



REF. 100122

**SUPER OFFER**  
**130 €**

### Specifications

Work pressure 4 to 5.5 bares
Particle size Up to 100 microns



## Gutta Percha cutter

The tip is flexible design to let you cut smoothly.

Press the power button and the tip can heat up to 500°C immediately. Release the button and the tip will cool down quickly.




REF. 100349

**90 €**



## X-ray viewers flat



**CHECK PRICES**

DIMENSIONS	VISION AREA	OUTDOOR AREA	USE	WEIGHT
Ref. 100131 DIN A4	300 × 210 mm	360 × 240 × 20 mm	For wall and tabletop	1.7 kg
Ref. 100132 SRA3	430 × 355 mm	510 × 495 × 50 mm	For wall and tabletop	3.1 kg
Ref. 100134 DIN A2	580 × 400 mm	700 × 520 × 33 mm	For wall	4 kg

## Silicone Mixer

- ✓ Suitable for 99% of silicones.



### Specifications

Height 290 mm	Width 170 mm
Depth 250 + 150 mm	Weight 7.1 kg
Power 150 W	Push force 450 kg ± 10
Environmental temp. 5-40°C	

REF. 100360

**930 €**

## Composite Heater

Device to keep up to 8 composite syringes warm.



REF. 080077

**SUPER OFFER**  
**225 €**

### Specifications

Height 135 mm	Width 97 mm
Depth 110 mm	Weight 0.7 kg
Temperature 40-50-70°C	Holes Ø 16 mm

## Portable micromotor brushless intramatic

Autonomy 10 hours



### CONTROL BOX

Height 45 mm	Width 105 mm
Depth 145 mm	Weight 0,4 kg

### REF. 100282

Height 20 mm	Width 20 mm	Depth 155 mm
Weight 0,2 kg	Power 30W	Torque 250 gf-cm
Speed 30000 rpm		

REF. 100282

**583 €**



# SPECIAL OFFER | Clinic

From 1 February to 31 May, 2025




## Magnifying glasses



Lightweight, handy and with excellent panoramic and clear vision. The magnifier is essential for clinical and laboratory procedures. Manufactured with CR-39 resin, higher resistance against abrasion and excellent and excellent transparency.

- ✓ Available in 2.5x or 3.5x magnification.
- ✓ Packing dimensions: (WxDxH) 25.5 x 16.5 x 3.5 cm.
- ✓ Weight: 65g.



NEW

### LAUNCHING OFFER

Augmentation 2,5 x

REF. 100076

36 €

Augmentation 3,5 x

REF. 100077

36 €



## Safety glasses

Vent valve

REF. 100014

10 €



Anti-fog

REF. 100011

15 €

UV rays

REF. 100013

15 €

Great Vision

REF. 100012

15 €

## Face protection

Folding + 10 units



REF. 100022

SUPER OFFER  
12 €

## Face protection

Frame + 10 units



REF. 100021

SUPER OFFER  
13 €

## Finishing box Eclipse



- ✓ Designed and manufactured in Spain.
- ✓ Coupling for jeb type chipper.
- ✓ Ergonomic working position.
- ✓ LED lamps.
- ✓ Easy cleaning.

REF. 080216


378 €



275 x 300 x 360 | 2,8 kg

## Flask

for duplicating with alginate



- ✓ The inner measurements can be reduced by using rubber pieces



REF. 020206

97 €

## Vibrador redondo

Ø120 mm



- ✓ Made of white colour anodized aluminium, and rests on suction pads.
- ✓ Recess of 30 mm for inserting the tray.
- ✓ It is equipped with a potentiometer, a selector switch (+,-) for stepped adjustment of vibrations, and an "ON" light.



Specifications

Height 100 mm  
Surfaces Ø 120 mm  
Weight 2.5 kg

REF. 080030

277 €

## Polymerizing pneumatic Pot




Inox

Specifications

Pressure máx. 2 bar  
Temp. máx. 60°C  
Capacity 6 l

REF. 030420

300 €



235 mm  
420 mm  
245 mm  
Inner dimensions: Ø 220x170 mm height  
2.4 kg

Double safety valve

## Water bath for wax plates

For intraoral records

- ✓ Stainless tank and lid.

Specifications

Height 155 mm  
Width 235 mm  
Depth 205 mm  
Weight 1.7 kg  
Temperature 30-75°C  
Capacity 1 l



REF. 080078

402 €

# SPECIAL OFFER | Clinic

From 1 February to 31 May, 2025



MESTRA

## Trolley

- ✓ Made of stainless steel.
- ✓ Pull-out tray (495 x 345 mm)



NEW

REF. 100451

**SUPER OFFER**  
**145 €**

885-1340 mm  
550 mm  
495 mm  
4,5 kg

Inox



## Trolley

- ✓ Made of stainless steel.
- ✓ Pull-out tray (495 x 345 mm)



NEW

REF. 100452

**SUPER OFFER**  
**145 €**


785-1220 mm  
550 mm  
495 mm  
4 kg

Inox



## Scanner trolley

NEW




875 mm  
465 mm (tray)  
550 mm (wheels)  
600 mm  
TONS -- kg



REF. 100450

**380 €**



## Electric mobile trolley

- ✓ Made of steel sheet and enameled white.
- ✓ Two sockets for equipment connection.
- ✓ Three trays.
- ✓ Four wheels (smooth and balanced movement).

760 mm  
290 mm  
400 mm  
12,5 kg  
TONS

REF. 100305

**297 €**



## Electric mobile trolley

- ✓ Made of steel sheet and enameled in white.
- ✓ Three trays.
- ✓ Three sockets to power the equipment.
- ✓ Four wheels with smooth and balanced movement.
- ✓ Weight supported per tray: 5 kg

REF. 100307

**SUPER OFFER**  
**216 €**



760 mm  
300 mm  
460 mm  
11 kg  
TONS



## Electric mobile trolley

- ✓ Made of steel sheet and enameled in white.
- ✓ Three trays.
- ✓ Three sockets, one drawer.
- ✓ Four wheels with smooth and balanced movement.
- ✓ Weight supported per tray: 5 kg.

REF. 100308

**SUPER OFFER**  
**366 €**



765 mm  
340 mm  
400 mm  
15 kg  
TONS



## Electric mobile trolley

- ✓ Made of steel sheet and enameled white.
- ✓ Two sockets for equipment connection.
- ✓ Four wheels (smooth and balanced movement).
- ✓ Three glass trays.

REF. 100306

**SUPER OFFER**  
**420 €**

740 mm  
360 mm  
460 mm  
20 kg  
TONS



## Working table

- ✓ Metallic table, wood countertop.
- ✓ Working chair  
Ref. 100300
- ✓ Polishing lathe and two trays.
- ✓ LED lamp.
- ✓ Micromotor.

Gift:

- ✓ Magnifying glass
- ✓ Magnetic bur holder

### Work chair

Allows four different regulations.  
Wheels with safety brakes  
(EN 1335-2:2009)



REF. 100301K2

**1.867 €**

REF. 100300

**266 €**





# SPECIAL OFFER | Clinic

From 1 February to 31 May, 2025



MESTRA

## Sprayer nebulizer

Using the appropriate products, large surfaces can be efficiently disinfected in just a few minutes (up to 100 m<sup>2</sup> in just 10 minutes).



### Specifications

Dimension	250 x 335 x 570 mm
Weight	2.7 kg
Power	1000 W
Capacity	5 l

REF. 100520

329 €

## Extraoral suction system

Greatly reduce the risk of cross contamination between professionals and patients and protect respiratory health.



### Specifications

Height (arm not included)	755 mm
Width	330 mm
Depth	320 mm
Weight (arm included)	30 kg
Power	1000 W
Arm length	1.3 m

REF. 100347

1.118 €

## UV cabinet sanitizer

The effectiveness of ultra violet rays is known to eliminate germs from any surface.  
Ideal for clinics and laboratories.



### Specifications

Dimension	300 x 290 x 200 mm
Weight	4.5 kg
Power	4 W

REF. 100502

186 €

## UV cabinet sanitizer

The effectiveness of ultra violet rays is known to eliminate germs from any surface.  
Ideal for small objects such as mobile phones, masks, key rings, coins, jewelry, etc.

- ✓ Highly recommended for personal use.
- ✓ With drying function.



### Specifications

Dimension	120 x 296 x 195 mm
Weight	1 kg
Power	150 W

REF. 100500

53 €

## UV cabinet sanitizer

The inside of the cabinet is made entirely of easy-to-clean stainless steel.  
It contains 10 levels of trays in which instruments can be stored.  
The deposited elements will be subjected to the action of UV rays.



### Specifications

Dimension	500 x 350 x 260 mm
Weight	8 kg
Power	15 W

REF. 100504

239 €

## Air purifier (60 - 70 m<sup>2</sup>)

It acts on the stale air in any room, eliminating dust, dirt and germs.



### Specifications

Dimension	570 x 370 x 265 mm
Weight	8.5 kg
Power	72 W

REF. 100514

406 €

## Air purifier (30 - 35 m<sup>2</sup>)

It acts on the stale air in any room, eliminating dust, dirt and germs.



### Specifications

Dimension	536 x 360 x 206 mm
Weight	6.3 kg
Power	45 W

REF. 100512

312 €

## Air purifier (25 m<sup>2</sup>)

It acts on the stale air in any room, eliminating dust, dirt and germs.



### Specifications

Dimension	303 x 329 x 170 mm
Weight	4 kg
Power	36 W

REF. 100510

179 €



Semiadjustable arcon type articulator + Case

### Arquímedes Professional A7 Plus

REF. 010099

200 €

Robust structure made of aluminum, ensuring greater stability. Design that provides greater visibility.

- ✓ Stabilization system of the condilar guide movement with silicone connection.
- ✓ Support for upper branch in open position
- ✓ Fixed intercondylar distance (110 mm)
- ✓ Curve condylar guide
- ✓ Regulable condylar guide
- ✓ Regulable Bennet angle
- ✓ Central lock



Semiadjustable arcon type articulator + Case

### Arquímedes Protrusal A7 Plus E

REF. 010103

SUPER OFFER  
285 €

It has exclusive aluminum condylar guides with millimeter adjustment.

- ✓ Stabilization system of the condilar guide movement with spring.
- ✓ Support for upper branch in open position.
- ✓ Fixed intercondylar distance (110 mm).
- ✓ Regulable condylar guide.
- ✓ Regulable Bennet angle.
- ✓ Flat condylar guide.
- ✓ Bilateral locking system.




Semiadjustable arcon type articulator + Case

### 4000 S

REF. 010107

240 €

- ✓ Adjustable intercondylar distance (expanding spindle).
- ✓ Curve condylar guide.
- ✓ Regulable condylar guide.
- ✓ Regulable Bennet angle.
- ✓ Stabilization system of the condilar guide movement with O-ring.



Arcon type articulator

### Arquímedes Junior A7 Fix

REF. 010098

112 €

- ✓ Fixed intercondylar distance (110 mm).
- ✓ Condylar guide set at average 30°.
- ✓ Bennet angle set at average 15°.
- ✓ Curve condylar guide.
- ✓ Magnetic stabilization system of the condilar guide movement (built-in).
- ✓ Support for upper branch in open position.



## Standard face bow

REF. 010102-50

95 €

The standard face bow is an apparatus used to record the spatial position of the patient's upper dental arch in relation to the skull and transfer the record to the articulator.

Transfer of the patient record is performed by setting the complete face bow in the upper frame of the articulator.



CE Medical device

## Elite face bow

REF. 010100-55

198 €

The new elite face bow brings an articulated and fixation system for the fork that allows a faster registration process for the professional and comfortable for the patient.

It has a jig transfer assembly that eliminates the needs for face bow in the assembling process of the casts in the articulator, turning the process faster, more practical and stable.

Its nazium relator with adjustable height, including an infra-orbital indicator, turns the registration process more precise.



CE Medical device

Is equipped with a sliding system which allows smooth lateral movement (opening/closing) and facilitates the overall use.

The ear pieces are anatomical, removable and chemically sterilized.





# SPECIAL OFFER | Clinic

From 1 February to 31 May, 2025



MESTRA

## Standard Professional Set

- ✓ Articulator A7 PLUS
- ✓ Standard face bow
- ✓ Case

REF. 010102

**SUPER OFFER**  
**265 €**



## Elite Professional Set

- ✓ Articulator A7 PLUS
- ✓ Élite face bow
- ✓ Case

REF. 010101

**379 €**



## Standard Protrusional Set

- ✓ Articulator A7 PLUS-E
- ✓ Standard face bow
- ✓ Case

REF. 010104

**SUPER OFFER**  
**330 €**



## Elite Protrusional Set

- ✓ Articulator A7 PLUS-E
- ✓ Élite face bow
- ✓ Case

REF. 010106

**SUPER OFFER**  
**450 €**



## Standard 4000 S Set

- ✓ Articulator 4000S
- ✓ Standard face bow
- ✓ Case

REF. 010107S

**331 €**



## Elite 4000 S Set

- ✓ Articulator 4000S
- ✓ Élite face bow
- ✓ Case

REF. 010107E


**418 €**




# Dry piston compressors

Family of compressors consisting of two sizes of 80 and 160 litres/minute to which a dryer with automatic drain can be coupled. Extremely quiet operation (60 dB), easy maintenance (dry piston) and excellent reliability. They are supplied with two pressure gauges (boiler and service pressure), outlet pressure regulator and complete set of valves for draining. The service pressure can reach 8 bars.


**Not suitable for milling machines.**



CCertificate by TÜV SÜD  
TYPE CERTIFICATE  
N° 12202200491258001  
Directive 2014/29/UE  
(Pressure vessels)



1. Ref. 110340 is also available with HEPA filter (H14) Ref. 110340H
2. The activated carbon filter Ref. 110300-53 is available as an accessory for all compressors.



REF. 110340

160 l/min  
(with dryer)

**1.974 €**




REF. 110310

80 l/min  
(with dryer)

**1.334 €**

REFERENCE	110310	110340
Height	650 mm	750 mm
Width	390 mm	680 mm
Length	540 mm	480 mm
Weight	37 kg	66 kg
Power	750 W	1500 W
Pressure	5-8 bar	5-8 bar
Boiler capacity	30 l	50 l
N° heads	1	2

## LED lamps and lamp base




**Best seller!**

REF. 100288

**SUPER OFFER**  
**90 €**

- ✓ 84 large LEDs
- ✓ Color temperature 6000 K
- ✓ Tactile dimmable light intensity.
- ✓ Specially rigid arm.
- ✓ With clamp for fastening to the table.



REF. 100290

**SUPER OFFER**  
**90 €**

- ✓ 117 LEDs
- ✓ Color temperature 5600-6000 K
- ✓ With clamp for fastening to the table.



**NEW**

REF. 100289

**SUPER OFFER**  
**90 €**

- ✓ 90 large LEDs
- ✓ Color temperature 6000-7000 K
- ✓ Tactile dimmable light intensity.
- ✓ With clamp for fastening to the table.



**Magnifying glass**

3 diopter lens

REF. 100287

5 diopter lens

REF. 100287-5x

**SUPER OFFER**  
**76 €**

- ✓ 90 white LEDs
- ✓ Ø 190 mm lens
- ✓ Clamp to fix it to the table




**NEW**


REF. 100289V

**108 €**

- ✓ 140 LEDs.
- ✓ Colour temperature 3000/6500 K
- ✓ 3 types of light: cold, neutral and warm.
- ✓ Tactile dimmable light intensity.
- ✓ With clamp for fastening to the table.



Swivels up to 180°



## Accessories for lamp and magnifying glass

REF. 100287 - REF. 100288



**NEW**

REF. 100288-10  
Wheeled stand

**57 €**



REF. 100288-02

Clamp support



REF. 100288-03

Wall bracket



REF. 100288-04

Screwable  
table support

25 cm



REF. 100288-07

Desktop stand



**MESTRA**

Talleres Mestraitua SL  
(+34) 944 530 388  
info@mestra.es - www.mestra.es  
48150 SONDIKA (Spain)

Offer valid unless typographical error or out of stock. The manufacturer reserves the right to make changes without prior notice.



Management  
System  
ISO 9001:2015

www.tuv.com  
ID 1100018003