

Promotions

OCTOBER – DECEMBER 2026



ENDODONTICS

CanalPro CalciSeal **NEW!**

p. 11



RESTORATION

BRILLIANT Crios

p. 14



ROTARY INSTRUMENTS

DIATECH Diamond burs

p. 30

 **COLTENE GROUP**

PRODUCT RANGES

COLTENE



RESTORATION

BRILLIANT range



ENDODONTICS

CanalPro range
HyFlex range



ROTARY INSTRUMENTS

DIATECH range
KENDA range



PROSTHETICS

AFFINIS
PRESIDENT The Original
Jet Blue Bite



TREATMENT AUXILIARIES

INFECTION CONTROL

ROEKO range
SciCan range



CanalPro X-Move

Full Freedom of Movement

The wireless CanalPro X-Move endomotor provides two movement modes, a full rotating and a reciprocation motion for all common file systems. The motor stands for simple handling and great flexibility.



Contra angle:
360 degrees
rotation



Torque range:
0.4 – 5.0 Ncm



Speed range:
100 – 2500 rpm



Adjustable
reciprocating angle:
10° interval



Adjustment range
20°–400°



Mini insulated
contra angle with
diameter of 8mm
and a height of
9.7mm



CanalPro X-Move

65002739

SPECIAL OFFER

12 FREE PRODUCTS

for the purchase of a **CanalPro X-Move** motor:

HyFlex EDM OGSF Sequence 25mm, 4 pcs. (60036142)	×5
HyFlex EDM Finishing Set 25mm, 3 pcs. (60019364)	×1
HySolate Latex Dental Dam Black medium, 36 pcs. (60031352)	×1
ROEKO GuttaFlow bioseal Standard Set 5ml (60019560)	×1
HyFlex EDM OGSF Gutta points 30/.04, 60 pcs. (60036504)	×1
HyFlex EDM Gutta points Assortment, 60 pcs. (60022705)	×1
HyFlex EDM OGSF Paper points 30/.04, 100 pcs. (60036503)	×1
HyFlex EDM Paper points Assortment, 100 pcs. (60022744)	×1

CanalPro X-Move pack + FREE products

SALESBUNDLEXM1



For more
information, scan
this QR code

HyFlex – OGSF Sequence



PREDICTABLE RESULTS

by applying the OGSF standard protocol




FAST PREPARATION OF DAILY CASES

thanks to the optimized match of tip and taper along the sequence



HIGHER PATIENT SAFETY

as a result of high material fracture resistance and following the OGSF standard protocol



— Working length

CanalPro

Irrigation and rinsing

CanalPro SyringeFill Holder

60031348 **163.21 € + VAT**

1 × Bottle holder for wall mounting
3 × CanalPro SyringeFill Cap

CanalPro SyringeFill Station

60031347 **269.44 € + VAT**

1 × Bottle holder incl. base
3 × CanalPro SyringeFill Cap



-25%
discount on
CanalPro solutions

CanalPro NaOCl 3 %

1 × 500ml

65002866 ~~40.90 € + VAT~~ **30.68 € + VAT**



CanalPro NaOCl 6 %

1 × 500ml

65002867 ~~40.90 € + VAT~~ **30.68 € + VAT**

CanalPro SyringeFill Cap

60019662 **11.68 € + VAT**

1 × SyringeFill Cap



CanalPro EDTA 17 %

1 × 500ml

60019651 ~~85.22 € + VAT~~ **63.92 € + VAT**




CanalPro Slotted end Tip

60011167 **42.11 € + VAT**

100 × Slotted-end tips with sideways opening Calibre 27 - 0.4 mm

60011168 **50.68 € + VAT**

100 × Slotted-end tips with sideways opening Calibre 30 - 0.3 mm



CanalPro Color Syringes

15.03 € + VAT

50 × Color Syringes 5ml

18.17 € + VAT

50 × Color Syringes 10ml

	5ml	10ml
●	60019321	60011173
●	60019322	60011174
●	60019323	60011175
●	60019324	60011176

CanalPro CHX 2 %

1 × 500ml

60019658 ~~65.62 € + VAT~~ **49.22 € + VAT**



CanalPro Side Port Tips

60011169 **107.54 € + VAT**

100 × Side-Port tips. Calibre 27 - 0.4 mm

60011170 **117.32 € + VAT**


100 × Side-Port tips. Calibre 30 - 0.3 mm



ENDODONTIC RETREATMENT

HyFlex Remover

The essential piece
for your retreatment



ONE FILE

to remove the root filling material



RESPECTS

the root canal anatomy



EFFICIENT

without any solvent



SAFE

thanks to non-cutting tip



30/.07

Maximum torque: **2.5 N.cm**

Rotation speed: **400-800 rpm**

Heat treatment: **C.Wire**



Obturation removal
up to WL -3 mm

HyFlex Remover

5 × B5 HyFlex Remover n30/7% Classics

67.97 € + VAT

19 mm	60023649
23 mm	60023648

NEW SHADE AVAILABLE

REINFORCED CAD/CAM COMPOSITE FOR PERMANENT RESTORATION

BRILLIANT Crios

BRILLIANT Crios is the ideal choice for single-tooth restorations, both in the anterior and the posterior region. This includes all conventional indications such as **inlays**, **onlays**, **crowns** and **veneers**.



INLAYS



ONLAYS



CROWNS



VENEERS



A-SILICONE IMPRESSION MATERIALS

AFFINIS PRESIDENT The Original

NOVEMBER
-20%
discount on AFFINIS
and PRESIDENT
The Original



For more
information on
the AFFINIS range,
scan this QR code



For more
information on the
PRESIDENT range,
scan this QR code

PUTTY

AFFINIS putty

115.08 € + VAT

2 × 300ml Base and Catalyst
2 × Dosage spoon



PRESIDENT The Original putty

93.07 € + VAT

2 × 300ml Base and Catalyst
2 × Dosage spoon



PUTTY	AFFINIS	PRESIDENT The Original
putty	-	 60019942
putty soft	 6530	 60019943
putty soft fast	 6531	 60022166
putty super soft	 6535	 60019944

RESORBABLE HAEMOSTATIC GELATINE SPONGES
FOR RELIABLE AND EFFICIENT BLEEDING CONTROL

roeko Gelatamp



	Gelatamp with colloidal silver			Gelatamp white			Gelatamp forte		
Strengths	Prolonged antibacterial activity due to colloidal silver			Promotes hemostasis and stabilizes the clot			Immediate absorption thanks to its special sponge-like structure		
Composition	Animal gelatin + colloidal silver			Animal gelatin			Animal gelatin		
Colour	Brown			White			White		
Dimensions	14 × 7 × 7 mm			10 × 10 × 10 mm			10 × 10 × 4 mm		
Absorption	 ABSORPTION  COLLOIDAL SILVER > 35 times its weight			 ABSORPTION 50 times its weight			 ABSORPTION 60 times its weight		
Resorption	4 weeks			4 weeks			4 weeks		
Use	Dry			Dry or moistened with sterile saline solution			Dry		
Box of 20 pcs.		274002	22.46€ + VAT	Box of 24 pcs.		60052842	18.82€ + VAT	Box of 24 pcs.	
Box of 50 pcs.		274007	51.67€ + VAT						
Box of 10 double blisters (20 pcs.)		274008	24.96€ + VAT						

HySolate Dental Dam latex



Isolated operating
field and optimal
vision



Safety and
protection

-20%
discount on a
selection of items



HySolate Latex Dental Dam - Powder free, low protein

Size	Colour	Thickness								Quantity				
		Thin		Medium		Heavy		X-heavy						
152 × 152 mm	 Green	60034418	~~21.90 € + VAT~~	17.52 € + VAT	60034419	~~24.32 € + VAT~~	19.46 € + VAT	60034420	~~27.90 € + VAT~~	22.31 € + VAT	60034421	~~28.75 € + VAT~~	23 € + VAT	×36
	 Black Edition				60031352	~~27.99 € + VAT~~	22.39 € + VAT	60031353	~~32.09 € + VAT~~	25.67 € + VAT	60031354	~~33.08 € + VAT~~	26.46 € + VAT	
	 Blue				60034442	~~24.32 € + VAT~~	19.46 € + VAT	60034443	~~27.90 € + VAT~~	22.31 € + VAT				
	 Light	60034427	~~21.90 € + VAT~~	17.52 € + VAT	60034428	~~24.32 € + VAT~~	19.46 € + VAT	60034429	~~27.90 € + VAT~~	22.31 € + VAT	60034430	~~28.75 € + VAT~~	23 € + VAT	
	 Fiesta (assortment)	60034433		21.90 € + VAT	60034436		25.26 € + VAT	60034438		27.99 € + VAT				
127 × 127 mm	 Green	60034414		21.90 € + VAT	60034415		24.32 € + VAT	60034416		27.90 € + VAT	60034417		28.75 € + VAT	×52
	 Blue				60034441		24.32 € + VAT							
	 Light	60034424		21.90 € + VAT	60034425		24.32 € + VAT	60034426		27.90 € + VAT				
	 Fiesta (assortment)	60034434		21.90 € + VAT	60034435		25.26 € + VAT	60034437		27.99 € + VAT				

HySolate Dental Dam non-latex

-20%
discount on HySolate
SyntX Dam

NEW!



HySolate
Dental Dam
Non-Latex



HySolate
Flexi Dam
Non-Latex



HySolate
Flexi Dam
Non-Latex
Framed



HySolate
SyntX Dam
Non-Latex

HySolate Dental Dam Non-Latex

Size	Colour	Medium	Quantity
127 x 127 mm	Green	H09928	42.13 € + VAT ×15
152 x 152 mm	Green	H09105	42.13 € + VAT ×15
	Green	H09106	218.24 € + VAT ×75

HySolate Flexi Dam Non-Latex

Size	Colour	Ref.	Quantity
152 x 152 mm	Blue	H09947	48.02 € + VAT ×30
	Green	H09946	48.02 € + VAT ×30
	Purple	H09945	48.02 € + VAT ×30

HySolate Flexi Dam Non-Latex Framed

Size	Colour	Medium	Quantity
100 x 105 mm	Purple framed	H00750	83.06 € + VAT ×20

HySolate SyntX Dam Non-Latex

Size	Colour	Medium	Quantity
127 x 127 mm	Blue	60037566	42.84 € + VAT ×30
	Blue	60037567	46.02 € + VAT ×30
152 x 152 mm	Blue	60037568	112.10 € + VAT ×75



For more
information, scan
this QR code

SciCan BRAVO G4

Chamber autoclave available
in 3 sizes: **17L** **22L** **28L**



**G4
technology**



**Smart & efficient
features**



**Various filling
options**

SciCan BRAVO G4
Dim. 60 × 48 × 46 cm

17 l

22 l

28 l

Ask your
COLTENE
distributor for
a price quote.



Cycle	Cycle Type	17 L		22 L		28 L	
		Cycle Time	Max Load	Cycle Time	Max Load	Cycle Time	Max Load
134°C Universal	B	42 min	6kg	46 min	7.5kg	56 min	9kg
134°C Hollow Unwrapped	S	35 min	6kg	39 min	7.5kg	44 min	9kg
134°C Solid Wrapped	S	33 min	6kg	39 min	7.5kg	45 min	9kg
134°C Prion	B	56 min	6kg	60 min	7.5kg	70 min	9kg
121°C Universal	B	58 min	6kg	63 min	7.5kg	69 min	9kg
134°C Helix Test	n.a.	20 min	n.a.	24 min	n.a.	24 min	n.a.
Vacuum Test	n.a.	18 min	n.a.	18 min	n.a.	18 min	n.a.
Custom (134°C or 121°C)	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.



For more
information, scan
this QR code