

get it!
NSK Special Offers



AUTUMN/WINTER - 31 AUGUST - 31 DECEMBER 2026

NEW



Oral Hygiene

Varios Combi Pro2

- An evolution in ultrasonic and powder therapy
- Removable Powder Chambers
- No routine servicing required
- Wireless multi-function foot-control
- Superbly comfortable water temperature, for both ultrasonic & powder



Minimally Invasive Therapy

Create a healthier future,
with you and M.I.T

Oral Surgery



VarioSurg4

Restorative and Prosthetics



Implantology



Surgic Pro2

Contra-angles & Air Turbines

Contra-angle Ti-Max Z99L

The ultimate in smooth cutting.
Designed to be durable.

THE ONLY HANDPIECE EVER TO
EARN A PERFECT SCORE!

5.0
RATING



Dental Product Shopper Evaluation

optic 1:5 **1.235€***

Ti-Max Z99L

~~1.487€*~~

REF. CA13140001

- Head size: Ø8.9 x H 13.4 mm
- Solid titanium body with DURAGRIP® coating
- Ceramic ball bearings
- Microfilter
- Two-way spray function
- Anti-Heat System
- 100° Angled head



DURAPower CHUCK®
reducing the risk of bur slippage



Quattro Spray



100° Angled head
for easy access



Optimised gear structure – for smooth cutting

Contra-angle S-Max M94L

Faster and smoother
zirconia removal.

optic 1:4 **906€***

S-Max M94L

~~949€*~~

REF. CA13170001

- Maximum speed 160,000 min
- Head size: Ø 9.9 x H 14.5 mm
- Stainless Steel
- For FG Burs (Ø1.6)
- Cellular Glass Optics
- Ceramic ball bearings



Exceptional Durability

Durability has improved by 2.7 times, delivering reliable performance when cutting hard materials such as zirconia, with excellent wear resistance.

Ti-Max Z Turbines

- Solid titanium body with DURAGRIP® coating • Cellular glass optics • Ceramic ball bearing • Clean Head System
- Quattro spray • Quick Stop • Simple cartridge change • Z990/Z890: Microfilter, DURAPower CHUCK® • 2 year warranty

Standard-head (Ø 12.1 x H 12.7 mm)

optic

Mini-head (Ø 10.6 x H 10.8 mm)

optic

Micro-head (Ø 9.0 x H 10.8 mm)

optic

Ti-Max Z990L

REF. PA23730001

- for NSK PTL coupling

Ti-Max Z990KL

REF. PA23750001

- for KaVo®MULTIflex® LUX coupling

Ti-Max Z990SL

REF. PA23770001

- for Sirona® quick coupling

Ti-Max Z990WL

REF. PA23870001

- for W&H® Roto Quick® coupling

1.051€*

~~1.427€*~~

Ti-Max Z890L

REF. PA23720001

- for NSK PTL coupling

Ti-Max Z890KL

REF. PA23740001

- for KaVo®MULTIflex® LUX coupling

Ti-Max Z890SL

REF. PA23760001

- for Sirona® quick coupling

Ti-Max Z890WL

REF. PA23860001

- for W&H® Roto Quick® coupling

1.051€*

~~1.427€*~~

Ti-Max Z micro

REF. PA2357

- for NSK PTL coupling

950€*

~~1.011€*~~

Ti-Max Z micro KL

REF. PA2358

- for KaVo®MULTIflex® LUX coupling

Ti-Max Z micro SL

REF. PA2359

- for Sirona® quick coupling

Ti-Max Z micro WL

REF. PA23610001

- for W&H® Roto Quick® coupling

1.025€*

~~1.138€*~~

NOTE: Ti-Max Z890/Z990 Air Turbines series are also available for Bien Air® couplings at the same promotion price

Contra-angles



Ti-Max Z

- Extremely small dimensions (e.g. Head size Ø 8.9 x H 13.4 mm) • Solid titanium body with scratch resistant DURAGRIP® coating
- Cellular glass optics • Integrated microfilter • Clean Head System • Push Button Chuck • 3 year warranty



optic 1:5

Ti-Max Z95L
REF. C1034

- Optic contra-angle
- 1:5 Increasing
- Quattro spray
- Ceramic ball bearings
- Two-way spray function

1.221€*
1.449€*



optic 1:1

Ti-Max Z25L
REF. C1038

- Optic contra-angle
- 1:1 Direct drive
- Two-way spray function

828€*
927€*



optic 4:1

Ti-Max Z15L
REF. C1039

- Optic contra-angle
- 4:1 Reduction
- Two-way spray function

945€*
1.073€*



optic 1:4,2

Ti-Max Z45L
REF. C1064

- Optic contra-angle
- 1:4.2 Increasing
- Two-way spray function

1.328€*

Ti-Max X

- Solid titanium body with scratch resistant DURAGRIP® coating • Cellular glass optics • Clean Head System • 2 year warranty



optic 1:5

Ti-Max X95L
REF. C600

- Optic contra-angle
- 1:5 Increasing
- Quattro spray
- Ceramic bearings

1.027€*
1.201€*



optic 1:1

Ti-Max X25L
REF. C601

- Optic contra-angle
- 1:1 Direct drive
- Single spray

754€*
885€*



optic 4:1

Ti-Max X15L
REF. C602

- Optic contra-angle
- 4:1 Reduction
- Single spray

863€*
995€*



optic 1:1

Ti-Max X65L
REF. H273

- Optic handpiece
- 1:1 Direct drive
- Single spray

896€*
976€*

Contra-angles & Air Turbines

S-Max M

• Stainless steel body • Cellular glass optics • Clean Head System • 2 year warranty




optic 1:5

S-Max M95L

REF. C1023

- Optic contra-angle
- 1:5 Increasing
- Quattro spray
- Ceramic ball bearing

906€*
949€*



optic 1:1

S-Max M25L

REF. C1024

- Optic contra-angle
- 1:1 Direct drive
- Single spray

645€*
714€*



optic 4:1

S-Max M15L

REF. C1025

- Optic contra-angle
- 4:1 Reduction
- Single spray

721€*
848€*



non optic 1:1

S-Max M65

REF. H1008

- Non-optic handpiece
- 1:1 Direct drive
- Single spray

490€*
520€*

• Stainless steel body • Cellular glass optics • Clean Head System • Quattro spray • 2 year warranty

Standard-head (Ø 12.1 × H 13.3 mm) optic



S-Max M900L

REF. P1254

- for NSK PTL coupling

688€*
733€*

S-Max M900KL

REF. P1258

- for KaVo®MULTIflex® LUX coupling

754€*
823€*

S-Max M900SL

REF. P1262

- for Sirona® quick coupling

754€*
823€*

S-Max M900WL

REF. P1266

- for W&H® Roto Quick® coupling

754€*
823€*

Mini-head (Ø 10.6 × H 12.4 mm) optic

S-Max M800L

REF. P1253

- for NSK PTL coupling

688€*
733€

S-Max M800KL

REF. P1257

- for KaVo®MULTIflex® LUX coupling

754€*
823€*

S-Max M800SL

REF. P1261

- for Sirona® quick coupling

754€*
823€*

S-Max M800WL

REF. P1265

- for W&H® Roto Quick® coupling

754€*
823€*


NOTE: S-Max M air turbines are also available for Bien Air® couplings at the same promotion price.

Ti-Max S970 Airscaler

- Titanium airscaler handpiece with DURACOAT® coating
- 3-level Power Control Ring
- Stable torque and constant power
- Circular optic illumination
- Wide selection of tips for perio, endo and scaling as well as restoration
- Frequency of oscillation: 5.800–6.200 Hz



Smooth and Easy 3-level
Power Control Ring



Unrestricted visibility



optic

Ti-Max S970L

REF. T1005

- for NSK PTL coupling

1.138€*

optic

Ti-Max S970KL

REF. T1013

- for KaVo®MULTIflex® LUX coupling

1.261€*

optic

Ti-Max S970SL

REF. T1015

- for Sirona® quick coupling

1.261€*



32
reasons to
buy the new
Varios Combi Pro2

NEW



Varios Combi Pro2 Basic Set

REF. Y1500505

6.129€*

Varios Combi Pro2 Basic Set contents:

- Varios Combi Pro2 Control Unit • VA2-LUX-HP Ti
- VA2-LUX-SC Ti • Prophy-Chamber2 • Prophy-Hose2
- Powder Handpiece (Prophy) • Wireless Foot-Control
- Hanger • Air Hose • Water Hose • Irrigation Bottle Set
- CLM Bottle Set • x1 FLASH pearl bottle
- Other accessories

Optional items:

Varios Combi Pro2
with water bottle kit

JET Nozzle
for Perio handpiece
REF. TA10770106
116€*

Perio Nozzle
for Perio handpiece
REF. T1077012
116€*

Prophy Nozzle
for Prophy handpiece
REF. T1074012
106€*

b.clean
Ready-to-use solution for disinfection and biofilm
prevention in water systems
1l bottles x 6
REF. BCL600
101€*



VALUE PACK P6 light

Varios Combi Pro2 Basic Set
REF. Y1500505

+ Varios Combi Pro2 Perio Set
Subgingival use
REF. Y1500470

6.480€*

~~7.414€*~~

SAVE 934€



VALUE PACK P6+C

Varios Combi Pro2 Basic Set
REF. Y1500505

+ Varios Combi Pro2 Perio Set
Subgingival use
REF. Y1500470

+ iCart Prophy2
Portable cart
REF. S9098

7.200€*

~~8.206€*~~

SAVE 1.006€



VALUE PACK P6+

Varios Combi Pro2 Basic Set
REF. Y1500505

+ Varios Combi Pro2 Perio Set
Subgingival use
REF. Y1500470

+ Varios Combi Pro2 Prophy Set
Supragingival use
REF. Y1500469

+ Titanium Ultrasonic-
Handpiece VA2-LUX
REF. EA1203X050

7.720€*

~~9.610€*~~

SAVE 1.890€



VALUE PACK P6+CW

Varios Combi Pro2 Basic Set
REF. Y1500505

+ Varios Combi Pro2 Perio Set
Subgingival use
REF. Y1500470

+ iCart Prophy2
Portable cart
REF. S9098

+ Water-Bottle-Kit
for iCart Prophy2
REF. Y1500719

7.820€*

~~8.850€*~~


SAVE 1.030€



Titanium slim
scaler handpiece



CLM mode for safe
reprocessing



Jet Nozzle for
Perio Handpiece



Water heating
system



Wireless multi-function
foot-control




Minimally Invasive Therapy

What is M.I.T.?

M.I.T. aligns with the benchmark in periodontal treatment – being anchored in the EFP S3 Level Clinical Guidelines and built around a patient-first philosophy. This modern approach empowers clinicians to deliver precise, predictable, and evidence-based care for stage I–III periodontitis.

By embracing this structured, research-driven framework, clinicians can elevate their practice, enhance patient outcomes, and align with global standards in periodontal care. This framework empowers clinicians to deliver effective, evidence-based care through three key pillars:

- Thorough Patient Assessment
- Evidence-Based Treatment
- Patient Empowerment



Create a healthier
future, with you and
M.I.T.



Prophy-Mate neo

SUPRAGINGIVAL AIR POLISHER



VALUE PACK P1

- Prophy-Mate neo

VARIOUS CONNECTION TYPES, SEE OVERVIEW
- + 2nd Prophy-Mate neo handpiece

PMN-HP

REF. T797010
- + FLASH pearl 4 x 300 g

REF. Y900698



1.215€*

1.530€*

SAVE 315€



Available connection types for Prophy-Mate neo:

For NSK PTL coupling	REF. Y135029	1.026€*
For KaVo® MULTiflex® LUX coupling	REF. Y135030	1.075€*
For Sirona® quick coupling	REF. Y135031	1.075€*
For W&H® Roto Quick® coupling	REF. Y135032	1.075€*
For Bien-Air® Unifix® coupling	REF. Y135033	1.075€*

MEGAPACK P3

Perio-Mate + Prophy-Mate neo

VALUE PACK P1 & P2 IN COMBINATION

2.437€*

3.294€*

SAVE 857€



Perio-Mate

BIOFILM ERASER



VALUE PACK P2

- Biofilm Eraser

Perio-Mate

VARIOUS CONNECTION TYPES, SEE OVERVIEW
- + 2nd handpiece

PER-HP

REF. T1043010
- + Glycine-based Powder

Perio Mate Powder 2 x 160 g

REF. Y900938



1.388€*

1.764€*

SAVE 376€



Available connection types for Prophy-Mate neo:

For NSK PTL coupling	REF. Y1002654	1.201€*
For KaVo® MULTiflex® LUX coupling	REF. Y1002655	1.266€*
For Sirona® quick coupling	REF. Y1002656	1.266€*
For W&H® Roto Quick® coupling	REF. Y1002657	1.266€*
For Bien-Air® Unifix® coupling	REF. Y1002658	1.266€*

Oral Hygiene

Powder Solution

FOR LONG LASTING, BEAUTIFUL AND HEALTHY TEETH



Calcium carbonate based powder

FLASH pearl

REF. Y900698

- 4 x 300g bottles
- Average particle diameter: 54µm
- For supragingival use
- For stain removal

97€*



Glycine-based powder

SOFT pearl

REF. Y1500721

- 4 x 250g bottles
- Average particle diameter: 65 µm
- For supragingival use
- Ideal for moderate staining.
- Water soluble

130€*



Glycine-based powder

Perio Mate Powder

REF. Y900938

- 2 x 160g bottles
- Average particle diameter: 25 µm
- For supragingival and subgingival use

77€*

Varios 970

ULTRASONIC SCALER

Varios 970 Complete Set contains:

- Control unit
- Handpiece incl. cord
- Foot control
- 3 scaling tips (G4, G8, G16)
- 2 irrigation bottles

Ultrasonic Scaler **non optic**

Varios 970

REF. Y1004160

2.721€*

Ultrasonic Scaler **optic**

Varios 970 LUX

REF. Y1004161

3.102€*



VALUE PACK P7

Ultrasonic Scaler **optic**

Varios 970 Lux

REF. Y1004161



Air Polisher for supragingival use

Prophy-Mate neo

3.089€*


4.177€*

SAVE 1.088€



Available connection types for Prophy-Mate neo:

For NSK PTL coupling	REF. Y135029	1.026€*
For KaVo® MULTIflex® LUX coupling	REF. Y135030	1.074€*
For Sirona® quick coupling	REF. Y135031	1.074€*
For W&H® Roto Quick® coupling	REF. Y135032	1.074€*
For Bien-Air® Unifix® coupling	REF. Y135033	1.074€*



NSK offers a wide choice of ultrasonic scaler tips for clinical applications.



iProphy Complete set

REF. Y1500015

- Cordless prophy handpiece with iFX57 head (4:1)
- Alternatively a PR-F head (optional extra) can be used that connects to disposable Prophy heads.

999€*

1.346€*

SAVE 347€



non optic 16:1

Ti-Max X57

REF. C619

- Contra-angle
- 16:1 Reduction
- Ultra-small head for screw-in cups and brushes

505€*



non optic 4:1

FX57m

REF. C1054

- Non-optic stainless steel contra-angle
- 4:1 Reduction
- Ultra-small head for screw-in cups and brushes
- Paste Defense System
- 360° Rotation

304€*

Oral Surgery & Implantology

iSD900

CORDLESS PROSTHETIC SCREWDRIVER

- Time-saving through automatic fastening of abutment screws
- Universally applicable for all major implant systems
- Torque-Calibration-System (TCS)
- Individually adjustable torque between 10 and 40 Ncm
- Selectable speeds (15, 20, 25/min)
- Ergonomic On/Off switch lever

Prosthetic screwdriver

iSD900

REF. Y1001358

1.689€*

~~1.893€*~~


SAVE **204€**



Osseo 100 / Osseo 100+

OSSEOINTEGRATION MONITORING DEVICE

Osseo 100 / 100+ measures implant stability to support decision making on when to load an implant. Especially important when working with shorter treatment time or managing risk patients. The uncomplicated procedure that measures ISQ allows the implant loading period to be planned in advance. The reconstruction of crowns and bridges can be monitored to optimise timing to decrease the risk for failures. Measurements can be made without unnecessary impact since the equipment does not come into physical contact with the implant or abutment.



Wireless Bluetooth®
connection with
Surgic Pro2

Osseo 100+

REF. Y1004176

1.899€*

~~2.816€*~~

SAVE **917€**

Tip: Get the Osseo 100+ for only 1.699*€ when you choose a surgery package with SurgicPro2 or VarioSurg 4.



Osseo 100

REF. Y1004175

1.699€*


~~2.481€*~~

SAVE **782€**



Only 3 steps
are necessary.
Learn more in
the video.

Advantages of Osseo 100+



Reduce
treatment time



Control and
maximum
safety



Advanced data
management
in combination
with
Surgic Pro2

MulTipeg not included, sold separately.

MulTipeg™

40€*



Reusable MulTipeg™

- For all major implant systems ¹⁾
- Tissue friendly, durable Titanium
- 20x autoclavable
- Optimal platform fit

¹⁾ There are different MulTipegs available especially designed for diverse implant systems and models.



Find your
MulTipeg

Oral Surgery & Implantology

Surgic Pro2

SURGICAL MICROMOTOR SYSTEM

Surgic Pro2 Complete Set contains:

- Control Unit
- SGL80M optic or SG80M non-optic micromotor
- X-SG20L optic contra-angle (20:1 Reduction) or SG20 non-optic contra-angle
- Wireless foot-control
- Irrigation tube (Pack of 3)
- Other accessories

Complete Set **non optic**

Surgic Pro2 NON-OPT

REF. Y1004196

3.937€*

Complete Set **optic**

Surgic Pro2 OPT

REF. Y1004195

5.166€*



VALUE PACK S6-B

Complete Set **non optic**
Surgic Pro2 NON-OPT
REF. Y1004196



+ 2nd surgical contra-angle (20:1)
S-Max SG20M **non optic**
REF. CA13020001



+ Surgical handpiece (1:1)
SGS-ES **non optic**
REF. H264



4.699€*

5.274€*

SAVE **575€**

VALUE PACK S7-B

Complete Set **optic**
Surgic Pro2 OPT
REF. Y1004195



+ 2nd surgical contra-angle (20:1)
Ti-Max X-SG20ML **optic**
REF. CA13030001



+ Surgical handpiece (1:1)
Ti-Max X-SG65L **optic**
REF. H1009



6.499€*

7.363€*

SAVE **864€**

Osseo 100+

REF. Y1004176



Buy one of the Surgic Pro2 packages S6-B or S7-B and you will receive the Osseo 100+ for only

1.699€*

2.816€*

SAVE **1.117€**



Connecting with Surgic Pro2 is very easy. Learn more in the video.

Surgical Handpieces

non optic **1:1**
SGS-ES
REF. H264

- 1:1 Direct drive
- Microsurgery-handpiece

504€*



optic **1:1**
Ti-Max X-SG65L
REF. H1009

1.000€*

non optic **1:1**
Ti-Max X-SG65
REF. H1038

621€*



NEW
non optic **20:1**
S-Max SG20M
REF. CA13020001

- 20:1 Reduction
- Stainless steel body
- Up to 80 Ncm Torque
- Mini-head

833€*



NEW
optic **20:1**
Ti-Max X-SG20ML
REF. CA13030001

- Optic contra-angle
- 20:1 Reduction
- Titanium body
- Up to 80 Ncm Torque
- Mini-head

1.197€*

Oral Surgery

VarioSurg 4

ULTRASONIC BONE SURGERY SYSTEM WITH LED



optic

VarioSurg 4

REF. Y1500706

6.886€*

optic

VarioSurg 4 w/o FC

REF. Y1500712

5.830€*

For combination with Surgic Pro2, without foot control.

VarioSurg 4 Complete Set:

- Control Unit
- LED Handpiece with 2m cord
- FC-86 Foot Control
- Sterilisation Cassette
- Handpiece Stand
- Irrigation tube (3 pcs.) and other accessories
- Basic H-S Kit (H-SG1, SG3, SG5, SG6D, SG7D, SG11 and Tip holder)

Main Features:

- 3 Modes (SURG, ENDO, PERIO)
- Selectable Burst Mode (3 frequencies)
- 9 Programms (SURG x 5, ENDO x 2, PERIO x 2)
- Adjustable LED Light Intensity
- Automatic Cleaning Mode
- Feedback-Function
- Auto-Tuning Function
- Dynamic Link Feature

Applications:

- Oral Surgery
- Bone surgery
- Scraper
- Sinus lifts
- Extraction
- Sinus Membrane Tool
- Implant preparation
- Perio
- Scaling
- Perio (Root Planing)
- Maintenance
- Retrograde Endo



User-friendly display

VALUE PACK V2

Complete Set with wireless foot control

VarioSurg 4

REF. Y1500706

The set without foot control is not available as part of the V2 value pack.



+ 3 additional Varios Tips of your choice

Example



+ Tip Holder

REF. Z221086

59€*

6.879€*

SAVE up to 1.136€



Powerful Handpiece with a Super-Slim Design



SG12

REF. ZT02010001

190€*



NSK Offers a wide choice of Vario Surg Surgical Tips for clinical applications.

Oral Surgery & Implantology

Surgic Pro2 × VarioSurg 4

TWO DEVICES. ONE SOLUTION. UNLIMITED POTENTIAL



Easy connection of Surgic Pro2 and VarioSurg 4

VarioSurg 4 features built-in Bluetooth® functionality, allowing it to be easily linked with the surgical micromotor system, Surgic Pro2. This enables the operation of both units with one single wireless foot control. Easily switch between the two systems, ensuring a smooth and efficient workflow.



VALUE PACK S10+

Complete Set

Surgic Pro2 OPT optic

REF. Y1004195



2nd optic contra-angle

+ Ti-Max X-SG20ML optic

REF. CA13030001



Surgical optic handpiece (1:1)

+ Ti-Max X-SG65L optic

REF. H1009



Ultrasonic-Bone-Surgery-System

+ VarioSurg 4 w/o FC optic

(without foot control)

REF. Y1500712



Mount to combine the two devices

+ Link Stand 3

REF. ZA16230001



11.139€*

13.601€*

SAVE **2.462€**

VALUE PACK S11+

Complete Set

Surgic Pro2 OPT optic

REF. Y1004195



2nd optic contra-angle

+ Ti-Max X-SG20ML optic

REF. CA13030001



Surgical optic handpiece (1:1)

+ Ti-Max X-SG65L optic

REF. H1009



Ultrasonic Bone Surgery System

+ VarioSurg 4 w/o FC optic

(without foot control)

REF. Y1500712



Special cart

+ iCart Duo

REF. S9090



12.199€*

14.158€*

SAVE **1.959€**



Osseo 100+

REF. Y1004176

Buy one of the Surgic Pro2/VarioSurg 4 value packs S10+ or S11+ and you will receive the Osseo 100+ (page 9) for only

1.699€*

2.816€*

SAVE **1.117€**



Connecting with Surgic Pro2 is very easy. Learn more in the video.

Oral Surgery & Implantology

Surgic Pro

SURGICAL MICROMOTOR FOR IMPLANTOLOGY, ORAL & MAXILLOFACIAL SURGERY

- Powerful torque (up to 80 Ncm)
- Small & lightweight micromotor
- LED illumination (32.000 LUX)
- Outstanding durability & reliability
- Micromotor autoclavable & thermo-disinfectable



Standard Sets

non optic

Surgic Pro NON-OPT

REF. Y1003587

Complete Set non-optic incl.
contra-angle (non-detachable)
SG20 (20:1 Reduction)

3.577€*

optic

Surgic Pro OPT

REF. Y1003586

Complete Set with LED incl.
optic contra-angle (non-detachable)
X-SG20L (20:1 Reduction)

4.771€*

optic

Surgic Pro+ OPT-D

REF. Y1003585

Complete Set with LED and data-storage
facility incl. detachable optic contra-angle
X-DSG20L (20:1 Reduction)

5.129€*

VALUE PACK S1-A

Complete Set

Surgic Pro NON-OPT

REF. Y1003587



2nd surgical contra-angle (20:1)

+ S-Max SG20 non optic

REF. C1010



Surgical handpiece (1:1)

+ SGS-ES non optic

REF. H264



4.259€*

4.914€*

SAVE **655€**

VALUE PACK S2-A

Complete Set

Surgic Pro OPT

REF. Y1003586



2nd surgical contra-angle (20:1)

+ Ti-Max X-SG20L optic

REF. C1003



Surgical handpiece (1:1)

+ Ti-Max X-SG65L optic

REF. H1009



5.949€*

6.968€*

SAVE **1.019€**

VALUE PACK S3-A

Complete Set

Surgic Pro+ OPT-D

REF. Y1003585



2nd surgical contra-angle (20:1)

+ Ti-Max X-DSG20L optic

REF. C1068



Surgical handpiece (1:1)

+ Ti-Max X-SG65L optic

REF. H1009



6.339€*

7.452€*

SAVE **1.113€**

Mobile Dentistry

Dentalone

THE MOBILE DENTAL UNIT

Dentalone

REF. S9070

Equipment:

- Brushless NLX nano micromotor with LED light
- Varios 170 LED multifunction ultrasonic scaler
- 3-way air/water/spray syringe with LED light
- Saliva suction
- Large illuminated control display
- Instrument carrier with automatic instrument selection
- Foot control
- Power supply required

13.516€*



VALUE PACK M1

Mobile dental unit

Dentalone

REF. S9070



Top Offer

Optic contra-angle (1:5)



S-Max M95L

optic



REF. C1023

9.890€*

14.465€*

SAVE **4.575€**

VIVAace

THE PORTABLE DENTAL UNIT

Whether in home care, in schools or on humanitarian missions, the mobile dental unit from NSK is very easy to transport and can be taken anywhere. (All you need is an electric supply.)

Weight: only 8.5 kg!



VIVAace

REF. Y1003773

VIVAace is equipped as follows:

- Brushless LED micromotor NLX plus with endo functions
- Multifunctional ultrasonic scaler Varios 170 with LED
- 3-Function syringe without light
- Saliva suction

11.673€*



VALUE PACK M2

Complete Set

VIVAace

REF. Y1003773



Top Offer

Optic contra-angle (1:5)



S-Max M95L

optic



REF. C1023

8.490€*

12.622€*

SAVE **4.132€**

VIVA MATE G5

PORTABLE MICROMOTOR WITH LIGHT FOR MOBILE USE

Wherever this mobile micromotor/control unit combination is used, in mobile dentistry or in the dental lab: thanks to its strong batteries it can be used totally independent even without a plug for up to 2 hours. With a rotation spectrum of 3,000 to 30,000 rpm it can be used for all types of dental applications when used with the right handpiece or contra-angle (not included in the set).

optic

VIVA MATE G5 Set

REF. Y1002801

1.820€*



PRESTO AQUA LUX

LUBE-FREE AIR TURBINE WITH WATER COOLING AND LED

- Speed: 320,000 rpm
- Cooling spray individually adjustable
- Minimal noise and vibration
- No lubrication required
- 360° free-swivel LED handpiece
- Unique dust-proof mechanism
- Integrated LED light

optic
PRESTO AQUA LUX
REF. Y1001151

- Lab turbine with LED

2.240€*
~~2.399€*~~

SAVE **159€**

non optic
PRESTO AQUA II
REF. Y150023

- Lab turbine without LED

1.735€*
~~1.836€*~~

SAVE **101€**



ULTIMATE XL

BRUSHLESS MICROMOTOR



Sets with Torque-handpiece
• (8,7 Ncm)

1.886€*

Sets with Compact-handpiece
• (6,0 Ncm)

1.580€*

- Speed range 1,000 – 50,000 rpm
- Lightweight, ergonomic handpiece
- Excellent, smooth running
- Patented dust protection mechanism
- Automatic speed control
- Auto-cruise function



Sets with Torque-handpiece

Knee control unit REF Y141495	Foot control unit REF Y141513	Table top control unit REF Y141471	Tower-type control unit REF Y141504
----------------------------------	----------------------------------	---------------------------------------	--

Sets with Compact-handpiece

Knee control unit REF Y141496	Foot control unit REF Y141514	Table top control unit REF Y141472	Tower-type control unit REF Y141505
----------------------------------	----------------------------------	---------------------------------------	--

VOLVERE i7

MICROMOTOR SYSTEM

VOLVERE i7 RM
REF. Y1002890

- Laboratory micromotor with lab handpiece

1.229€*



- Compact design • Automatic speed-control • Foot (pedal) or manual speed control
- Silent, smooth running • Speed range 1,000 to 35,000 rpm



VOLVERE i7 E
REF. Y1003523
• ISO-E type micromotor, without handpiece

1.101€*



VOLVERE i7 VR-EB
REF. Y1002891
• Typ ISO-E with lab handpiece

1.288€*

Contra-angles

FXm series

• Stainless steel body




non optic 4:1

FX15m

REF. C1053001

- 4:1 Reduction
- Max Speed: 10,000 min⁻¹
- For CA burs (ø2.35)
- Ball Bearing
- Push Button Chuck
- Clean Head System

377€*



non optic 1:1

FX25m

REF. C1052001

- 1:1 Direct Drive
- Max Speed: 40,000 min⁻¹
- For CA burs (ø2.35)
- Ball Bearing
- Push Button Chuck
- Clean Head System

261€*



non optic 4:1

FX57m

REF. C1054001

- 4:1 Reduction
- Max Speed : 5,000 min⁻¹
- 360° Rotation
- Paste Defense System
- External Water Spray
- For screw-in cups & brushes

304€*



non optic 1:1

FX65m

REF. H1014001

- 1:1 Direct Drive
- Max Speed: 40,000 min⁻¹
- External Water Spray
- For HP burs (ø2.35)
- Clean Head System

278€*



non optic 16:1

FX75m

REF. C1055001

- 16:1 Reduction
- Miniature Head
- Max Speed: 1,250 min⁻¹
- 360° Rotation
- For engine files
- Push Button Chuck

404€*



non optic 1:1

FX20m

REF. CA13190001

- 1:1 Direct Drive
- FG burs (ø2.35)
- Clean Head System
- Push Button chuck

253€*



non optic 4:1

FX50m

REF. HA12050001

- 4:1 Reduction
- For Disposable Prophy head (DPA)
- Max. Speed: 2.500 min⁻¹

325€*



non optic 4:1

FX55m

REF. CA13160001

- 4:1 Reduction
- For EVA-/Profin-/Orthofiles Tip
- Max Speed: 10,000 min⁻¹
- Vertical Reciprocating (Stroke Length 1.4 mm)

351€*



non optic 4:1

FX58m

REF. CA13210001

- 4:1 Reduction
- Max. Speed: 2.500 min⁻¹
- For screw-in cups & brushes
- Paste Defense System
- External Water Spray

304€*



non optic 16:1

FX70m

REF. CA13250001

- 16:1 Reduction
- Push Button chuck
- For Hand Files / 60°Twist
- Max. Speed: 2.500 min⁻¹

407€*



non optic 16:1

FX74m

REF. CA13200001

- 16:1 Reduction
- Push Button chuck
- For Engine Files (ø2.35) / 80° Twist

407€*

NEW



CONTACT INFORMATION FOR SOUTHEASTERN EUROPE



Robert Jurcevic

GENERAL MANAGER SALES AND MARKETING
CENTRAL & EASTERN EUROPE

✉ robert@nsk-nakanishi.co.jp

☎ +385 91 600 8092

* All prices are excluding VAT. All prices valid until 31st of December 2026.
Changes, errors, misprints etc. reserved. Contact your local NSK dealer for NSK's attractive offers.

NSK Europe GmbH

www.nsk-dental.com

Elly-Beinhorn-Str. 8 - 65760 Eschborn - Germany - Tel: +385 91 600 8092 - e-mail: robert@nsk-nakanishi.co.jp