



get it!

NSK Special Offers

AUTUMN / WINTER · 01. SEPTEMBER – 31. DECEMBER 2025

NEW



Oral Hygiene

Varios Combi Pro2

Oral Surgery



VarioSurg4

**Restorative and
Prosthetics**



Implantology



Surgic Pro2

Contra-angles & Air Turbines

Contra-angle

Ti-Max Z99L **NEW**

The ultimate in smooth cutting.
Designed to be durable.


optic 1:5 **1.199€***

Ti-Max Z99L 1.444€*
REF. CA13140001

- Head size: Ø8.9 x H 13.4 mm
- Solid titanium body with DURAGRIP® coating
- Ceramic ball bearings
- Microfilter
- Two-way spray function
- Anti-Heat-System
- 100° Angled head



DURAPOWER CHUCK®
reducing the risk of bur slippage



Quattro Spray



100° Angled head
for easy access



Optimized gear structure – for smooth cutting

Contra-angle

S-Max M94L **NEW**

Faster and smoother
zirconia removal

optic 1:4 **880€***

S-Max M94L 921€*
REF. CA13170001

- Maximum speed 160,000 min
- Head size: Ø 9.9 x H 14.5 mm
- Stainless Steel
- For FG Burs (Ø1.6)
- Cellular Glass Optics
- Ceramic ball bearings
- Clean Head System
- Quattro spray



Exceptional Durability

The durability has improved by 2.7 times, making it well-suited for zirconia's high hardness while also offering excellent wear resistance.

Ti-Max Z Air Turbines

• Solid titanium body with DURAGRIP® coating • Cellular glass optics • Ceramic ball bearing • Clean Head System • Quattro spray • Quick Stop • Simple cartridge change • Z990/Z890: Microfilter, DURAPOWER CHUCK® • 2 years warranty

Standard-head (Ø 12,1 x H 12,7 mm)

optic



Ti-Max Z990L

REF. PA23730001

- for NSK PTL coupling

Ti-Max Z990KL

REF. PA23750001

- for KaVo®MULTiflex® LUX coupling

Ti-Max Z990SL

REF. PA23770001

- for Sirona® quick-coupling

Ti-Max Z990WL

REF. PA23870001

- for W&H® Roto Quick® coupling

1.020€*

1.385€*

Mini-head (Ø 10,6 x H 10,8 mm)

optic



Ti-Max Z890L

REF. PA23720001

- for NSK PTL coupling

Ti-Max Z890KL

REF. PA23740001

- for KaVo®MULTiflex® LUX coupling

Ti-Max Z890SL

REF. PA23760001

- for Sirona® quick-coupling

Ti-Max Z890WL

REF. PA23860001

- for W&H® Roto Quick® coupling

1.020€*

1.385€*

Micro-head (Ø 9,0 x H 10,8 mm)

optic



Ti-Max Z micro

REF. PA2357

- for NSK PTL coupling

922€*

982€*

Ti-Max Z micro KL

REF. PA2358

- for KaVo®MULTiflex® LUX coupling

Ti-Max Z micro SL

REF. PA2359

- for Sirona® quick-coupling

Ti-Max Z micro WL

REF. PA23610001

- for W&H® Roto Quick® coupling

997€*

1.105€*

NOTE: Ti-Max Z890/Z990 Air Turbines series are also available for Bien Air® couplings at the same promotion price.

Contra-angles



Ti-Max Z

• Extremely small dimensions (e.g. Head size Ø 8,9 x H 13,4 mm) • Solid titanium body with scratch resistant DURAGRIP® coating • Cellular glass optics • Integrated microfilter • Clean Head System • Ultra-push chuck • 3 years warranty



optic 1:5

Ti-Max Z95L
REF. C1034

- Optic contra-angle
- 1:5 Increasing
- Quattro spray
- Ceramic ball bearings
- Two-way spray function

1.185€*
1.407€*



optic 1:1

Ti-Max Z25L
REF. C1038

- Optic contra-angle
- 1:1 Direct drive
- Two-way spray function

804€*
900€*



optic 4:1

Ti-Max Z15L
REF. C1039

- Optic contra-angle
- 4:1 Reduction
- Two-way spray function

917€*
1.042€*



optic 1:4,2

Ti-Max Z45L
REF. C1064

- Optic contra-angle
- 4:1 Reduction
- Two-way spray function

1.289€*

Ti-Max X

• Solid titanium body with scratch resistant DURAGRIP® coating • Cellular glass optics • Clean Head System • 2 years warranty



optic 1:5

Ti-Max X95L
REF. C600

- Optic contra-angle
- 1:5 Increasing
- Quattro spray
- Ceramic ball bearings

997€*
1.166€*



optic 1:1

Ti-Max X25L
REF. C601

- Optic contra-angle
- 1:1 Direct drive
- Single spray

732€*
859€*



optic 4:1

Ti-Max X15L
REF. C602

- Optic contra-angle
- 4:1 Reduction
- Single spray

838€*
966€



optic 1:1

Ti-Max X65L
REF. H273

- Optic handpiece
- 1:1 Direct drive
- Single spray

870€*
948€*

Contra-angles & Air Turbines

S-Max M

• Stainless steel body • Cellular glass optics • Clean Head System • 2 years warranty



optic 1:5

S-Max M95L

REF. C1023

- Optic contra-angle
- 1:5 Increasing
- Quattro spray
- Ceramic ball bearing

880€*
921€*



optic 1:1

S-Max M25L

REF. C1024

- Optic contra-angle
- 1:1 Direct drive
- Single spray

626€*
693€*



optic 4:1

S-Max M15L

REF. C1025

- Optic contra-angle
- 4:1 Reduction
- Single spray

700€*
823€*



non optic 1:1

S-Max M65

REF. H1008

- Non-optic handpiece
- 1:1 Direct drive
- Single spray

476€*
505€*

• Stainless steel body • Cellular glass optics • Clean Head System • Quattro-spray • 2 years warranty

Standard head (Ø 12,1 × H 13,3 mm) optic



S-Max M900L

REF. P1254

- for NSK PTL coupling

668€*
712€*

S-Max M900KL

REF. P1258

- for KaVo®MULTiflex® LUX coupling

732€*
799€*

S-Max M900SL

REF. P1262

- for Sirona® quick-coupling

732€*
799€*

S-Max M900WL

REF. P1266

- for W&H® Roto Quick® coupling

732€*
799€*

Mini head (Ø 10,6 × H 12,4 mm) optic

S-Max M800L

REF. P1253

- for NSK PTL coupling

668€*
712€*

S-Max M800KL

REF. P1257

- for KaVo®MULTiflex® LUX coupling

732€*
799€*

S-Max M800SL

REF. P1261

- for Sirona® quick-coupling

732€*
799€*

S-Max M800WL

REF. P1265

- for W&H® Roto Quick® coupling

732€*
799€*


NOTE: S-Max M air turbines are also available for Bien Air® couplings at the same promotion price.

Ti-Max S970 Aircaler

- Titanium aircaler handpiece with DURACOAT® coating
- 3-level Power Control Ring
- Stable torque and constant power
- Circular optic illumination
- Wide selection of tips for perio, endo and scaling as well as restoration
- Frequency of oscillation: 5.800–6.200 Hz



Smooth and Easy 3-level Power Control Ring



Unrestricted visibility



optic

Ti-Max S970L

REF. T1005

- for NSK PTL coupling

1.105€*

optic

Ti-Max S970KL

REF. T1013

- for KaVo®MULTiflex® LUX coupling

1.224€*


optic

Ti-Max S970SL

REF. T1015

- for Sirona® quick-coupling

1.224€*



32 good reasons to buy the new Varios Combi Pro2



Titanium slim
scaler handpiece



CLM-mode for safe
reprocessing



Jet Nozzle for
Perio Handpiece



Water heating
system



Wireless multi-
function foot-control

VALUE PACK P6

Varios Combi Pro2 Basic Set

+ Varios Combi Pro2 Perio Set
Subgingival use

6.290€*

~~7.198€*~~

SAVE **908€**

VALUE PACK P6+C

Varios Combi Pro2 Basic Set

+ Varios Combi Pro2 Perio Set
Subgingival use

+ iCart Prophy2
Portable cart

6.990€*

~~7.967€*~~

SAVE **977€**

VALUE PACK P6+

Varios Combi Pro2 Basic Set

+ Varios Combi Pro2 Perio Set
Subgingival use

+ Varios Combi Pro2 Prophy Set
Supragingival use

+ Titanium Ultrasonic-
Handpiece
VA2-LUX-HP-Ti

7.490€*

~~9.330€*~~

SAVE **1.840€**

NEW



Varios Combi Pro2 Basic Set

REF. Y1500505

5.950€*

Varios Combi Pro2 Basic Set contains:

Varios Combi Pro2 Control Unit • VA2-LUX-HP Ti •
VA2-LUX-SC Ti • Prophy-Chamber2 • Prophy-Hose2 •
Powder Handpiece (Prophy) • Wireless Foot-Control •
Hanger • Air Hose • Water Hose • Irrigation Bottle Set •
CLM Bottle Set • 1 Bottle FLASH pearl • and
other accessories

VALUE PACK P6+CW

Varios Combi Pro2 Basic Set

+ Varios Combi Pro2 Perio Set
Subgingival use

+ iCart Prophy2
Portable cart

+ Watter-Bottle-Kit
for iCart Prophy2

7.590€*

~~8.592€*~~

SAVE **1.002€**



Learn more about
the new
ultrasonic & powder
therapy unit
Varios Combi Pro2:



Oral Hygiene

Prophy-Mate neo

SUPRAGINGIVAL AIR POLISHER

VALUE PACK P1

Air Polisher

Prophy-Mate neo

VARIOUS CONNECTION TYPES, SEE OVERVIEW



2nd Prophy-Mate neo handpiece

+ PMN-HP
REF. T797010



Cleaning powder, based on calcium

+ FLASH pearl (4 x 300g)
REF. Y900698



1.178€*

~~1.487€*~~

SAVE **309€**



Available connection types for Prophy-Mate neo::

For NSK PTL coupling	REF. Y135029	996€*
For KaVo® MULTiflex® LUX coupling	REF. Y135030	1.044€*
For Sirona® quick-coupling	REF. Y135031	1.044€*
For W&H® Roto Quick® coupling	REF. Y135032	1.044€*
For Bien-Air® Unifix® coupling	REF. Y135033	1.044€*

MEGAPACK P3

Perio-Mate + Prophy-Mate neo

VALUE PACK P1 & P2 IN COMBINATION

2.366€* ~~3.201€*~~ SAVE **835€**



Perio-Mate

BIOFILM ERASER

VALUE PACK P2

Biofilm Eraser

Perio-Mate

ARIOUS CONNECTION TYPES, SEE OVERVIEW



2nd handpiece

+ PER-HP
REF. T1043010



Glycin based Powder

+ Perio Mate Powder (2 x 160g)
REF. Y900938



1.348€*

~~1.714€*~~

SAVE **366€**



Available connection types for Perio-Mate:

For NSK PTL coupling	REF. Y1002654	1.166€*
For KaVo® MULTiflex® LUX coupling	REF. Y1002655	1.229€*
For Sirona® quick-coupling	REF. Y1002656	1.229€*
For W&H® Roto Quick® coupling	REF. Y1002657	1.229€*
For Bien-Air® Unifix® Kupplung	REF. Y1002658	1.229€*

Oral Hygiene

Powder Solution

FOR LONG LASTING, BEAUTIFUL AND HEALTHY TEETH



Calcium carbonate based powder

FLASH pearl

REF. Y900698

- 4 x 300g bottles
- Average particle diameter: 54µm
- For supragingival use
- For stain removal

94€*



Glycine based powder

SOFT pearl

REF. Y1500721

- 4 x 250g bottles
- Average particle diameter: 65 µm
- For supragingival use
- ideal for moderate staining.
- Water soluble

126€*



Glycine based powder

Perio Mate Powder

REF. Y900938

- 2 x 160g bottles
- Average particle diameter: 25 µm
- For supragingival and subgingival use

75€*

Varios 970

ULTRASONIC SCALER

- iPiezo engine for stable power output
- Extremely slim & lightweight handpiece for all indications
- Extensive tip range for prophylaxis, periodontics, endodontics and restorative treatments
- LED integrated
- 2 irrigation bottles + external water

Varios 970 Complete Set includes:

- Control unit
- Handpiece incl. cord
- Foot control
- 3 scaling tips (G4, G8, G16)
- 2 irrigation bottles

Ultrasonic Scaler **non optic**

Varios 970

REF. Y1004160

2.642€*

Ultrasonic Scaler **optic**

Varios 970 LUX

REF. Y1004161

3.012€*



NSK offers a wide choice of ultrasonic scaler tips for clinical applications.



VALUE PACK P7

Ultrasonic Scaler

Varios 970 Lux

REF. Y1004161

optic



Air Polisher for supragingival use

+ Prophy-Mate neo

VARIOUS CONNECTION TYPES, SEE OVERVIEW

2.999€*

4.056€*

SAVE 1.057€



Available connection types for Prophy-Mate neo:

For NSK PTL coupling	REF. Y135029	996€*
For KaVo® MULTIflex® LUX coupling	REF. Y135030	1.044€*
For Sirona® quick-coupling	REF. Y135031	1.044€*
For W&H® Roto Quick® coupling	REF. Y135032	1.044€*
For Bien-Air® Unifix® coupling	REF. Y135033	1.044€*



iProphy Complete set

REF. Y1500015

- Cordless prophy handpiece
- with iFX57 head (4:1)
- Alternatively a PR-F head (optional extra) can be used that connects to disposable Prophy heads.

999€*

SAVE 308€

Prophy Contra-angles

non optic 16:1

Ti-Max X57

REF. C619

- Contra-angle
- 16:1 Reduction
- Ultra-small head for screw-in cups and brushes

490€*



non optic 4:1

FX57m

REF. C1054

- Non-optic stainless steel contra-angle
- 4:1 Reduction
- Ultra-small head for screw-in cups and brushes
- Paste Defense System
- 360° Rotation

295€*



Oral Surgery & Implantology

iSD900

CORDLESS PROSTHETIC SCREW DRIVER

- Time-saving through automatic fastening of abutment screws
- Universally applicable for all major implant systems
- Torque-Calibration-System (TCS)
- Individually adjustable torque between 10 and 40 Ncm
- Selectable speeds (15, 20, 25/min)
- Ergonomic On/Off switch lever

Prosthetic screw driver

iSD900

REF. Y1001358

1.645€*

1.838€*

SAVE **193€**



Osseo 100 / Osseo 100+

OSSEOINTEGRATION MONITORING DEVICE

Osseo 100 measures implant stability to support decision making on when to load an implant. Especially important when working with shorter treatment time or managing risk patients. The uncomplicated procedure that measures ISQ allows the implant loading period to be planned in advance. The reconstruction of crowns and bridges can be monitored to optimize timing to decrease the risk for failures. Measurements can be made without unnecessary impact since the equipment does not come into physical contact with the implant or abutment.



Wireless Bluetooth®
connection with
Surgic Pro2

Osseo 100+

REF. Y1004176

1.895€*

2.734€*

SAVE **839€**

Tip: Get the Osseo 100+ for only €1.699* when you choose a surgery package with SurgicPro2 or VarioSurg4.



Only 3 steps
necessary.
Learn more in
the video.

Osseo 100

REF. Y1004175

1.695€*

2.380€*

SAVE **685€**

MulTipeg not included, sold separately.

MulTipeg™

39€*



Reusable MulTipeg™

For all major implant systems ¹⁾


- Tissue friendly, durable Titanium
- 20x autoclavable
- Optimal platform fit

¹⁾There are different MulTipees available especially designed for diverse implant systems and models.



Find your
MulTipeg

Advantages of Osseo 100+



Reduce
treatment time



Control
and
maximum
safety



Advanced data
management
in combination
with
Surgic Pro2

Oral Surgery & Implantology

Go beyond.

Surgic Pro2

SURGICAL MICROMOTOR SYSTEM

Surgic Pro2 Complete Set contains:

- Control Unit
- SGL80M optic or SG80M non-optic micromotor
- X-SG20L optic contra-angle (20:1 Reduction) or SG20 non-optic contra-angle
- Wireless foot-control
- Irrigation tube (Pack of 3)
- and other accessories

Complete Set **non optic**

Surgic Pro2 NON-OPT

REF. Y1004196

3.822€*

Complete Set **optic**

Surgic Pro2 OPT

REF. Y1004195

5.016€*



VALUE PACK S6-B

Complete Set **non optic**
Surgic Pro2 NON-OPT
REF. Y1004196



Additional surgical non-optic
contra-angle (20:1)

+ S-Max SG20 **non optic**
REF. C1010



Surgical handpiece (1:1)

+ SGS-ES **non optic**
REF. H264



4.562€*

5.120€*

SAVE **558€**

VALUE PACK S7-B

Complete Set **optic**
Surgic Pro2 OPT
REF. Y1004195



Additional surgical non-optic
contra-angle (20:1)

+ Ti-Max X-SG20L **optic**
REF. C1003



Surgical handpiece (1:1)

+ Ti-Max X-SG65L **optic**
REF. H1009



6.313€*

7.149€*

SAVE **836€**

Osseo 100+

REF. Y1004176

Buy one of the Surgic Pro2 packages S6-B
or S7-B and you will receive the Osseo 100+ for only

1.699€*

2.734€*

SAVE **1.035€**



Connecting
with Surgic Pro2
is very easy.
Learn more
in the video.

Oral Surgery & Implantology

VarioSurg 4

ULTRASONIC BONE SURGERY SYSTEM WITH LED



optic

VarioSurg 4

REF. Y1500706

6.685€*

optic

VarioSurg 4 w/o FC

REF. Y1500712

5.660€*

For combination with Surgic Pro2, without foot control.

VarioSurg 4 Complete Set:

- Control Unit
- LED Handpiece with 2m cord
- FC-86 Foot Control
- Sterilisation Cassette
- Handpiece Stand
- Irrigation tube (3 pcs.) and other accessoires
- Basic H-S Kit (H-SG1, SG3, SG5, SG6D, SG7D, SG11 and Tip holder)

Main Features:

- Major Power Boost
- 3 Modes (SURG, ENDO, PERIO)
- Selectable Burst Mode (3 frequencies)
- 9 Programms (SURG x 5, ENDO x 2, PERIO x 2)
- Adjustable LED Light Intensity
- Automatic Cleaning Mode
- Feedback-Function
- Auto-Tuning Function
- Dynamic Link Feature

Applications

- Oral Surgery
- Bone surgery
- Scraper
- Sinus lifts
- Extraction
- Sinus Membrane Tool
- Implant preparation
- Perio
- Scaling
- Perio (Root Planing)
- Maintenance
- Retrograde Endo



User-friendly display

VALUE PACK V2

Complete Set with wireless foot control

VarioSurg 4

REF. Y1500706

The set without foot control is not available as part of the V2 value pack.



+ 3 additional Varios Tips of your choice



Example

+ Tip Holder

REF. Z221086

57€*



6.685€*

SAVE up to 780 €



Powerful Handpiece with a Super-Slim Design



NSK Offers a wide choice of Vario Surg Surgical Tips for clinical applications.

NEW

SG12

REF. ZT02010001

184 €*



Oral Surgery & Implantology

Surgic Pro2 × VarioSurg 4

TWO DEVICES. ONE SOLUTION. UNLIMITED POTENTIAL



Easy connection of Surgic Pro2 and VarioSurg 4

VarioSurg 4 features built-in Bluetooth® functionality, allowing it to be easily linked with the surgical micromotor system, Surgic Pro2. This enables the operation of both units with one single wireless foot control. Easily switch between the two systems, ensuring a smooth and efficient workflow.



VALUE PACK S10+

Complete Set

Surgic Pro2 OPT optic

REF. Y1004195



2nd optic contra-angle

+ Ti-Max X-SG20L optic

REF. C1003



Surgical optic handpiece (1:1)

+ Ti-Max X-SG65L optic

REF. H1009



Ultrasonic-Bone-Surgery-System

+ VarioSurg 4 w/o FC optic

(without foot control)

REF. Y1500712



Mounting to combine the two devices

+ Link Stand 3

REF. ZA16230001



10.814€*

13.205€*

SAVE **2.391€**

VALUE PACK S11+

Complete Set

Surgic Pro2 OPT optic

REF. Y1004195



2nd optic contra-angle

+ Ti-Max X-SG20L optic

REF. C1003



Surgical optic handpiece (1:1)

+ Ti-Max X-SG65L optic

REF. H1009



Ultrasonic-Bone-Surgery-System

+ VarioSurg 4 w/o FC optic

(without foot control)

REF. Y1500712



Special cart

+ iCart Duo

REF. S9090



11.844€*

13.746€*

SAVE **1.902€**



Osseo 100+

REF. Y1004176

Buy one of the Surgic Pro2/VarioSurg4 value packs S10+ or S11+ and you will receive the Osseo 100+ (page 8) for only

1.699€*

2.734€*

SAVE **1.035€**



Connecting with Surgic Pro2 is very easy. Learn more in the video.

Oral Surgery & Implantology

Surgic Pro

SURGICAL MICROMOTOR FOR IMPLANTOLOGY, ORAL & MAXILLOFACIAL SURGERY

- Powerful torque (up to 80 Ncm)
- Small & lightweight micromotor
- LED illumination (32.000 LUX)
- Outstanding durability & reliability
- Micromotor autoclavable & thermo-disinfectable



Standard Sets

non optic

Surgic Pro NON-OPT

REF. Y1003587

Complete Set non-optic incl. contra-angle (non-dismountable) SG20 (20:1 Reduction)

3.473€*

optic

Surgic Pro OPT

REF. Y1003586

Complete Set with LED incl. optic contra-angle (non-dismountable) X-SG20L (20:1 Reduction)

4.632€*

optic

Surgic Pro+ OPT-D

REF. Y1003585

Complete Set with LED and data-storage facility incl. dismountable optic contra-angle X-DSG20L (20:1 Reduction)

4.980€*

VALUE PACK S1-A

Complete Set

Surgic Pro NON-OPT

REF. Y1003587



2nd surgical contra-angle (20:1)

+ S-Max SG20 non optic

REF. C1010



Surgical handpiece (1:1)

+ SGS-ES non optic

REF. H264



4.138€*

~~4.771€*~~

SAVE **633€**

VALUE PACK S2-A

Complete Set

Surgic Pro OPT

REF. Y1003586



2nd surgical contra-angle (20:1)

+ Ti-Max X-SG20L optic

REF. C1003



Surgical handpiece (1:1)

+ Ti-Max X-SG65L optic

REF. H1009



5.783€*

~~6.765€*~~

SAVE **982€**

VALUE PACK S3-A

Complete Set

Surgic Pro+ OPT-D

REF. Y1003585



2nd surgical contra-angle (20:1)

+ Ti-Max X-DSG20L optic

REF. C1068



Surgical handpiece (1:1)

+ Ti-Max X-SG65L optic

REF. H1009



6.154€*

~~7.235€*~~

SAVE **1.081€**

Surgical handpieces and contra-angles



optic **20:1**
Ti-Max X-DSG20L
REF. C1068

1.284€*

non optic **20:1**
Ti-Max X-DSG20
REF. C1067

1.044€*

- 20:1 Reduction
- Solid titanium body
- Max Torque 80 Ncm
- Dismantleable



optic **20:1**
Ti-Max X-SG20L
REF. C1003

1.162€*

- 20:1 Reduction
- Solid titanium body
- Max Torque 80 Ncm



non optic **20:1**
S-Max SG20
REF. C1010

809€*

- 20:1 Reduction
- Stainless steel body
- Max Torque 80 Ncm



non optic **20:1**
SGX-E20R
REF. HA1200

1.426€*

- For Osteotomy preparation, insertion of Zygoma implants for Maxillofacial, oral surgery
- 20:1 Reduction
- Stainless steel body
- Max Torque 80 Ncm



optic **1:1**
Ti-Max X-SG65L
REF. H1009

971€*

non optic **1:1**
Ti-Max X-SG65
REF. H1038

603€*

- 1:1 Direct Drive
- Solid titanium body



non optic **1:1**
SGS-ES
REF. H264

489€*

- 1:1 Direct Drive
- Micro Surgery Handpiece



non optic **1:1**
SGA-ES
REF. H263

664€*

- 1:1 Direct Drive
- 20° Angle Handpiece



optic **1:1**
Ti-Max X-SG25L
REF. C1011

799€*

- 1:1 Direct Drive
- Solid titanium body



optic **1:3**
Ti-Max Z-SG45L
REF. C1107

1.204€*

non optic **1:3**
Ti-Max Z-SG45
REF. C1108

939€*

- 1:3 Increasing
- Solid titanium body
- 45° Angle Head



optic **1:3**
Ti-Max X-SG93L
REF. C1004

1.193€*

non optic **1:3**
Ti-Max X-SG93
REF. C1007

950€*

- 1:3 Increasing
- Solid titanium body



non optic **20:1**
SGMS-ER20i
REF. Y110147

720€*

- 20:1 Reduction
- External and Internal Cooling
- Incl. 2 depth indicators
- Head and shank separable



non optic **20:1**
SGM-ER20i
REF. Y110127

609€*

- 20:1 Reduction
- External and Internal Cooling
- Head and shank separable

Mobile Dentistry

VIVA ace

THE PORTABLE DENTAL UNIT

Whether in home care, in schools or on humanitarian missions, the mobile dental unit from NSK is very easy to transport and can be taken anywhere. (All you need is an electric supply.)

VIVA ace

REF. Y1003773

VIVA ace is equipped as follows:

- Brushless LED micromotor NLX plus with endo functions
- Multifunctional ultrasonic scaler Varios 170 with LED
- 3-F syringe without light
- Saliva ejector

11.333€*

Weight: only 8.5 kg !



VALUE PACK M2

Complete Set
VIVA ace
REF. Y1003773



S-Max M95L
(Optic contra-angle, 1:5)

+ S-Max M95L optic



REF. C1023

11.342€*

12.254€*

SAVE 912€



VIVA MATE G5

PORTABLE MICROMOTOR WITH LIGHT FOR MOBILE USE

Wherever this mobile micromotor/control unit combination is used, in mobile dentistry or in the dental lab: thanks to its strong batteries it can be used totally independent even without a plug for up to 2 hours. With a rotation spectrum of 3,000 to 30,000 rpm it can be used for all types of dental applications when used with the right handpiece or contra-angle (not included in the set).

optic

VIVA MATE G5 Set

REF. Y1002801

1.767€*



PRESTO AQUA LUX

LUBE-FREE AIR-TURBINE WITH WATER COOLING AND LED

- Speed: 320,000 rpm
- Cooling spray individually adjustable
- Minimal noise and vibration
- No lubrication required
- 360° free-swivel LED handpiece
- Unique dust-proof mechanism
- Integrated LED-light

optic

PRESTO AQUA LUX

REF. Y1001151

- Lab turbine with LED

2.175€*
2.329€*

SAVE **154€**

non optic

PRESTO AQUA II

REF. Y150023

- Lab turbine without LED

1.688€*
1.783€*

SAVE **95€**



ULTIMATE XL

BRUSHLESS MICROMOTOR



Sets with Torque-handpiece

- (8,7 Ncm)

1.831€*

Sets with Compact-handpiece

- (6,0 Ncm)

1.534€*

- Speed range 1.000 – 50.000 rpm
- Lightweight, ergonomic handpiece
- Excellent, smooth running
- Patented dust protection mechanism
- Automatic speed control
- Auto-cruise function



Sets with Torque-handpiece

Knee control unit
REF Y141495

Foot control unit
REF Y141513

Table top
control unit
REF Y141471

Tower-type
control unit
REF Y141504

Sets with Compact-handpiece

Knee control unit
REF Y141496

Foot control unit
REF Y141514

Table top
control unit
REF Y141472

Tower-type
control unit
REF Y141505

VOLVERE i7

MICROMOTOR SYSTEM

- Compact design • Automatic speed-control • Foot (pedal) or manual speed control • Silent, smooth running • Speed range 1,000 to 35,000 rpm

VOLVERE i7 RM

REF. Y1002890

- Laboratory micromotor with lab handpiece

1.193€*



VOLVERE i7 E

REF. Y1003523

- ISO-E type micromotor, without handpiece

1.069€*



VOLVERE i7 VR-EB

REF. Y1002891

- Typ ISO-E with lab handpiece

1.250€*

FX Series

• Stainless steel body

non optic 4:1  FX15m REF. C1053001 • 4:1 Reduction • Max Speed: 10,000 min⁻¹ • For CA burs (ø2.35) • Ball Bearing • Push Button Chuck • Clean Head System 366€*	non optic 1:1  FX25m REF. C1052001 • 1:1 Direct Drive • Max Speed: 40,000 min⁻¹ • For CA burs (ø2.35) • Ball Bearing • Push Button Chuck • Clean Head System 253€*	non optic 4:1  FX57m REF. C1054001 • 4:1 Reduction • Max Speed : 5,000 min⁻¹ • 360° Rotation • Paste Defense System • External Water Spray • For screw-in cups & brushes 295€*	non optic 1:1  FX65m REF. H1014001 • 1:1 Direct Drive • Max Speed: 40,000 min⁻¹ • External Water Spray • For HP burs (ø2.35), disposable prophylaxis angle • Clean Head System 270€*	non optic 16:1  FX75m REF. C1055001 • 16:1 Reduction • Miniature Head • Max Speed: 1,250 min⁻¹ • 360° Rotation • For engine files • Push Button Chuck 392€*
--	--	--	---	---

non optic 1:1  FX20m REF. CA13190001 • 1:1 Direct Drive • FG burs • Clean Head System • Push-Button chuck 246€*	non optic 4:1  FX50m REF. HA12050001 • 4:1 reduction • Disposable Prophylaxis Angles (DPA) 316€*	non optic 4:1  FX55m REF. CA13160001 • 4:1 reduction • EVA-/Profin-/Orthofiles Tip • 1.4mm stroke length 341€*	NEW Available in September 
non optic 4:1  FX58m REF. CA13210001 • 4:1 reduction • Snap-on cups & brushes • Paste defence system 295€*	non optic 16:1  FX70m REF. CA13250001 • 16:1 reduction • Push-Button chuck • 60° left/right rotation; for Endo files with hand shaft (hand files) 395€*	non optic 16:1  FX74m REF. CA13200001 • 16:1 reduction • Push-Button chuck • 80° left/right rotation; for Endo files with contra-angle shaft (2.35mm) 395€*	

* All prices are excluding VAT. All prices valid until December 31st, 2025. Changes, errors, misprints etc. reserved. Contact your local NSK dealer for NSK's attractive offers.